





		Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	FANNIN UK LTD	Supplier	SUPERMAX HEALTHCARE LIMITED
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Aurelia Zero	Brand	Aurelia Quest - M	Brand	Dermagrip Ultra LS	Brand	Aurelia Blush - M
		MPC	92777	MPC	92897	MPC	D1102 - 24 - SIZE M	MPC	78887
		NPC	FTE2474	NPC	FTE2055	NPC	FTE2493	NPC	FTE2070
		Newton	6N	Newton	6N	Newton	No info available	Newton	6N
		EN374 tested	No info available	EN374 tested	Yes	EN374 tested	No info available	EN374 tested	Yes
		Unit of Issue	2000	Unit of Issue	1000	Unit of Issue	200	Unit of Issue	2000
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (e.g. multiple gloves falling out)	✗		✗		✗		✗	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★☆☆ (1.75)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)	
	User able to maintain dexterity and grip when wet	★★★☆☆ (1.75)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)	
	User retains sensitivity of touch	★★★☆☆ (1.75)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)	
	Colour of glove does not impede clinical activity	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	





		Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	HEALTHLINE	Supplier	HEALTHLINE
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Aurelia Transform 100 - M	Brand	Aurelia Transform - M	Brand	Newton/Handsaf	Brand	Bodyguards
		MPC	9889A7	MPC	98897	MPC	Gn91/m	MPC	GL891/m
		NPC	FTE2060	NPC	FTE2065	NPC	FTE1164	NPC	FTE1187
		Newton	6N	Newton	6N	Newton	6N	Newton	6N
		EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	2000	Unit of Issue	2000	Unit of Issue	200	Unit of Issue	1000
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✓		✓		✓		✗	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)	
	User able to maintain dexterity and grip when wet	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (1.75)	
	User retains sensitivity of touch	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (1.5)	
	Colour of glove does not impede clinical activity	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	


		Supplier	ROBINSON HEALTHCARE LIMITED	Supplier	MEDICARE PRODUCTS LTD	Supplier	ROBINSON HEALTHCARE LIMITED	Supplier	SEMPERIT
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Nytraguard Bluple	Brand	Nitrex Ultra	Brand	Nytraguard Sky Blue	Brand	Sempermed Sempercure
		MPC	9547	MPC	GN02M	MPC	9527	MPC	826781635
		NPC	FTE1832	NPC	FTE1477	NPC	FTE1848	NPC	FTG178
		Newton	6N	Newton	No info available	Newton	6N	Newton	9N
		EN374 tested	Yes	EN374 tested	No info available	EN374 tested	Yes	EN374 tested	No info available
		Unit of Issue	200	Unit of Issue	200	Unit of Issue	1000	Unit of Issue	100
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✓		✗		✓		✓	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2.25)		★★★ (1.75)		★★★★ (2.5)		★★★★ (2.25)	
	User able to maintain dexterity and grip when wet	★★★★ (2.25)		★★★ (1.75)		★★★★ (2.5)		★★★★ (2.25)	
	User retains sensitivity of touch	★★★★ (2.25)		★★★★ (2)		★★★★ (2.5)		★★★★ (2.25)	
	Colour of glove does not impede clinical activity	★★★ (1.75)		★★★ (2)		★★★ (2)		★★★★ (2.25)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	





		Supplier	MEDICARE PRODUCTS LTD	Supplier	UNIGLOVES (UK) LIMITED	Supplier	SHERMOND(DIV BUNZL RET SUP)	Supplier	UNIGLOVES (UK) LIMITED
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Nitrex Extra Sensitive	Brand	Yellow Pearl	Brand	Premier Sensitive	Brand	Kooltouch
		MPC	GN01M	MPC	GP0113 (M)	MPC	P2741	MPC	GM0043 (M)
		NPC	FTE1746	NPC	FTE2343	NPC	FTE2358	NPC	FTE2298
		Newton	6N	Newton	6N	Newton	No info available	Newton	6N
		EN374 tested	No info available	EN374 tested	Yes	EN374 tested	No info available	EN374 tested	Yes
		Unit of Issue	200	Unit of Issue	1000	Unit of Issue	200	Unit of Issue	1000
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✗		✗		✗	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★☆ (2.25)		★★★★☆ (2.25)		★★★★☆ (2.25)		★★★★☆ (2)	
	User able to maintain dexterity and grip when wet	★★★★☆ (2.25)		★★★★☆ (2.25)		★★★★☆ (2.25)		★★★★☆ (2)	
	User retains sensitivity of touch	★★★★☆ (2.25)		★★★★☆ (2)		★★★★☆ (2)		★★★★☆ (2.25)	
	Colour of glove does not impede clinical activity	★★★★☆ (2)		★★★★☆ (2)		★★★★☆ (2)		★★★★☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	UNIGLOVES (UK) LIMITED	Supplier	UNIGLOVES (UK) LIMITED	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	HEALTHLINE
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Blue Pearl	Brand	Unicare Nitrile	Brand	Aurelia Bold - M	Brand	Newton/Handsafe
		MPC	GP0013 (M)	MPC	GU003 (M)	MPC	73997	MPC	Gn90/m
		NPC	FTE2318	NPC	FTE2278	NPC	FTE2036	NPC	FTG302
		Newton	6N	Newton	6N	Newton	6N	Newton	6N
		EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	1000	Unit of Issue	2000	Unit of Issue	1000	Unit of Issue	200
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (e.g. multiple gloves falling out)	✗		✗		✗		✗	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (2)	
	User able to maintain dexterity and grip when wet	★★★★ (2)		★★★★ (1.75)		★★★★ (1.75)		★★★★ (1.75)	
	User retains sensitivity of touch	★★★★ (2.25)		★★★★ (2.25)		★★★★ (2)		★★★★ (2)	
	Colour of glove does not impede clinical activity	★★★★ (2)		★★★★ (2)		★★★★ (1)		★★★★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	





		Supplier	SEMPERIT	Supplier	UNIGLOVES (UK) LIMITED	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	HEALTHLINE
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Sempermed	Brand	Pink Pearl	Brand	Supermax high newton nitrile	Brand	Handsafe
		MPC	826768535	MPC	GP0053 (M)	MPC	96897	MPC	GN99/m
		NPC	FTE1767	NPC	FTE2328	NPC	FTE2032	NPC	FTG026
		Newton	6N	Newton	6N	Newton	9N	Newton	9N
		EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	200	Unit of Issue	1000	Unit of Issue	1000	Unit of Issue	100
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✗		✗		✗	
	Able to don easily on dry hands	✓		✗		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2)		★★★★ (2.25)		★★★★ (2.25)		★★★★ (2)	
	User able to maintain dexterity and grip when wet	★★★★ (1.75)		★★★★ (2.25)		★★★★ (2.25)		★★★★ (2)	
	User retains sensitivity of touch	★★★★ (2.25)		★★★★ (2.25)		★★★★ (1.75)		★★★★ (2)	
	Colour of glove does not impede clinical activity	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	




		Supplier	UNIGLOVES (UK) LIMITED	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	HEALTHLINE	Supplier	SUPERMAX HEALTHCARE LIMITED
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Black Pearl	Brand	Aurelia Perform - M	Brand	Healthline	Brand	Aurelia Amazing - m
		MPC	GP0033 (M)	MPC	92397	MPC	Gn830/m	MPC	92887
		NPC	FTE2313	NPC	FTE2075	NPC	FTE1803	NPC	FTE2007
		Newton	6N	Newton	6N	Newton	No info available	Newton	6N
		EN374 tested	Yes	EN374 tested	Yes	EN374 tested	No info available	EN374 tested	Yes
		Unit of Issue	1000	Unit of Issue	2000	Unit of Issue	1000	Unit of Issue	300
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✓		✓		✓		✓	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★☆☆ (2)		★★★★☆ (2.25)		★★★☆☆ (2)		★★★★☆ (2.25)	
	User able to maintain dexterity and grip when wet	★★★☆☆ (2)		★★★★☆ (2.25)		★★★☆☆ (2)		★★★★☆ (2.25)	
	User retains sensitivity of touch	★★★☆☆ (2)		★★★★☆ (2.25)		★★★☆☆ (2)		★★★★☆ (2.25)	
	Colour of glove does not impede clinical activity	★★☆☆☆ (1)		★★★★☆ (2)		★★★☆☆ (2)		★★★★☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	HEALTHLINE	Supplier	HEALTHLINE	Supplier	DAVIS SCHOTTLANDER & DAVIS LTD	Supplier	DAVIS SCHOTTLANDER & DAVIS LTD
NITRILE NON-STERILE EXAMINATION GLOVES				Awaiting photo					
		Brand	Handsafe	Brand	Handsafe	Brand	Soft Touch Flexible Nitrile	Brand	Schottlander Flexible Nitrile
		MPC	Gn83/m	MPC	GN560	MPC	0153-M	MPC	877 - M
		NPC	FTG156	NPC	FTG274	NPC	FTG332	NPC	FTE775
		Newton	6N	Newton	6N	Newton	6N	Newton	6N
		EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	100	Unit of Issue	100	Unit of Issue	200	Unit of Issue	100
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✗		✓		✓	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2)		★★★★ (2.5)		★★★★ (2.25)		★★★★ (2.25)	
	User able to maintain dexterity and grip when wet	★★★★ (2)		★★★★ (2.25)		★★★★ (2.25)		★★★★ (2)	
	User retains sensitivity of touch	★★★★ (2)		★★★★ (2.5)		★★★★ (2.5)		★★★★ (2.25)	
	Colour of glove does not impede clinical activity	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	





		Supplier	SHERMOND(DIV BUNZL RET SUP)	Supplier	WILLIAMS MEDICAL SUPPLIES PLC	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	DAVIS SCHOTTLANDER & DAVIS LTD
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Performer	Brand	Williams medical supplies	Brand	Aurelia Quest 2.2 - M	Brand	Super soft flexible nitrile
		MPC	P2751	MPC	D4703M	MPC	92937	MPC	0260-M
		NPC	FTE1449	NPC	FTE2116	NPC	FTE2080	NPC	FTE1807
		Newton	6N	Newton	6N	Newton	6N	Newton	6N
		EN374 tested	Yes	EN374 tested	No info available	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	200	Unit of Issue	2000	Unit of Issue	2000	Unit of Issue	200
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✓		✓		✗		✓	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★☆ (1.75)		★★★★☆ (2.25)		★★★★☆ (2.25)		★★★★☆ (2.25)	
	User able to maintain dexterity and grip when wet	★★★★☆ (1.75)		★★★★☆ (2.25)		★★★★☆ (2)		★★★★☆ (2.25)	
	User retains sensitivity of touch	★★★★☆ (1.75)		★★★★☆ (2.25)		★★★★☆ (2.25)		★★★★☆ (2.25)	
	Colour of glove does not impede clinical activity	★★★★☆ (2)		★★★★☆ (2)		★★★★☆ (2)		★★★★☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	MEDICARE PRODUCTS LTD	Supplier	HEALTHLINE	Supplier	FANNIN UK LTD	Supplier	FANNIN UK LTD
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Nitrex Extra Length	Brand	Newton/Handsafe	Brand	Dermagrip	Brand	Dermagrip Ultra LT Soft
		MPC	EGN08M	MPC	GN92/M	MPC	D1502-19	MPC	D1102 - 92
		NPC	FTE1751	NPC	FTE1445	NPC	FTE1205	NPC	FTE1731
		Newton	6N	Newton	6N	Newton	6N	Newton	6N
		EN374 tested	No info available	EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	100	Unit of Issue	200	Unit of Issue	200	Unit of Issue	200
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✗		✗		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✓		✗		✗		✗	
	Able to don easily on dry hands	✓		✗		✗		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (1.75)	
	User able to maintain dexterity and grip when wet	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (1.75)	
	User retains sensitivity of touch	★★★☆☆ (1.75)		★★★☆☆ (2)		★★★★☆ (2.25)		★★★☆☆ (2)	
	Colour of glove does not impede clinical activity	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	FANNIN UK LTD	Supplier	MEDICARE PRODUCTS LTD	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	SUPERMAX HEALTHCARE LIMITED
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Dermagrip Ultra LT (LITE) Nitrile	Brand	Nitrex Sensitive	Brand	Aurelia Robust - M	Brand	Aurelia Sonic 200
		MPC	D1502-20	MPC	GN03M	MPC	93897	MPC	93777
		NPC	FTE1297	NPC	FTE1071	NPC	FTE2045	NPC	FTE2468
		Newton	6N	Newton	6N	Newton	6N	Newton	6N
		EN374 tested	Yes	EN374 tested	No info available	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	200	Unit of Issue	200	Unit of Issue	1000	Unit of Issue	200
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✓		✓		✓		✓	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (1.75)		★★★☆☆ (2)	
	User able to maintain dexterity and grip when wet	★★★☆☆ (1.75)		★★★☆☆ (2)		★★★☆☆ (1.75)		★★★☆☆ (2)	
	User retains sensitivity of touch	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (1.75)		★★★★☆ (2.25)	
	Colour of glove does not impede clinical activity	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (1.75)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	




		Supplier	UNIGLOVES (UK) LIMITED	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	SEMPERIT	Supplier	DAVIS SCHOTTLANDER & DAVIS LTD
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	White Pearl	Brand	Aurelia Robust Plus - M	Brand	Sempercare Skin 2	Brand	Simply Soft Flexible Nitrile
		MPC	GP0023 (M)	MPC	63887	MPC	826793845	MPC	0995-M
		NPC	FTE2338	NPC	FTE2040	NPC	FTG337	NPC	FTG431
		Newton	6N	Newton	6N	Newton	6N	Newton	6N
		EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	100	Unit of Issue	1000	Unit of Issue	200	Unit of Issue	200
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✓		✓		✗		✗	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (1.75)		★★★☆☆ (2)	
	User able to maintain dexterity and grip when wet	★★★☆☆ (2)		★★★☆☆ (1.75)		★★★☆☆ (2)		★★★☆☆ (2)	
	User retains sensitivity of touch	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (1.75)		★★★☆☆ (2)	
	Colour of glove does not impede clinical activity	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	




		Supplier	SHERMOND(DIV BUNZL RET SUP)	Supplier	SHERMOND(DIV BUNZL RET SUP)	Supplier	HEALTHLINE	Supplier	UNIGLOVES (UK) LIMITED
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Premier performer	Brand	Healthgard	Brand	Safedon	Brand	Green Pearl
		MPC	P2761	MPC	PSCR0152	MPC	GX894M	MPC	GP0043 (M)
		NPC	FTG374	NPC	FTG412	NPC	FTG385	NPC	FTE2323
		Newton	6N	Newton	6N	Newton	6N	Newton	6N
		EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	200	Unit of Issue	200	Unit of Issue	300	Unit of Issue	1000
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		N/A no box available	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		no box available	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		no box available	
	Product should fit in a dispenser	✓		✓		✓		no box available	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✗		✓		no box available	
	Able to don easily on dry hands	✓		✓		✓		no box available	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2)		★★★★ (2)		★★★★ (2.25)		★★★★ (2)	
	User able to maintain dexterity and grip when wet	★★★★ (2)		★★★★ (2)		★★★★ (2.25)		★★★★ (2)	
	User retains sensitivity of touch	★★★★ (2)		★★★★ (2)		★★★★ (2.25)		★★★★ (2)	
	Colour of glove does not impede clinical activity	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		no box available	

		Supplier	UNIGLOVES (UK) LIMITED	Supplier	Unigloves (Uk) Limited	Supplier	HEALTHLINE	Supplier	ROBINSON HEALTHCARE LIMITED
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Violet Pearl	Brand	Stronghold	Brand	Safedon	Brand	Nytraguard Tough Nitrile
		MPC	GP0073 (M)	MPC	GM0063	MPC	GX994M	MPC	9532
		NPC	FTE2333	NPC	FTE2548	NPC	FTG457	NPC	FTE1851
		Newton	6N	Newton	No info available	Newton	6N	Newton	6N
		EN374 tested	Yes	EN374 tested	No info available	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	1000	Unit of Issue	100	Unit of Issue	125	Unit of Issue	1000
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	N/A no box available		N/A Stronghold		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	no box available		no box available		✓		✓	
	Opening strong enough to resist tearing from repeated removal	no box available		no box available		✓		✓	
	Product should fit in a dispenser	no box available		no box available		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	no box available		no box available		✓		✗	
	Able to don easily on dry hands	no box available		no box available		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2)		★★★★ (2)		★★★★ (2.2)		★★★★ (2)	
	User able to maintain dexterity and grip when wet	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (1.7)	
	User retains sensitivity of touch	★★★★ (2)		★★★★ (1.5)		★★★★ (2.2)		★★★★ (1.7)	
	Colour of glove does not impede clinical activity	★★★★ (2.5)		★★★★ (2)		★★★★ (1.8)		★★★★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	no box available		no box available		✓		✓	

		Supplier	HEALTHLINE	Supplier	UNIGLOVES (UK) LIMITED	Supplier	MEDICARE PRODUCTS LTD	Supplier	HEALTHLINE
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Healthline	Brand	Zero	Brand	Medicare Nitrex	Brand	Newton/Handsafe
		MPC	Gn85/m	MPC	GM0053 M	MPC	Gn05m	MPC	GN892M
		NPC	FTE999	NPC	FTE2308	NPC	FTG222	NPC	FTE1232
		Newton	6N	Newton	6N	Newton	9N	Newton	6N
		EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	100	Unit of Issue	1000	Unit of Issue	100	Unit of Issue	200
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✗		✗		✗	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★☆ (2.25)		★★★★☆ (2.25)		★★★★☆ (2)		★★★★☆ (2.25)	
	User able to maintain dexterity and grip when wet	★★★★☆ (2.25)		★★★★☆ (2.25)		★★★★☆ (2)		★★★★☆ (2.25)	
	User retains sensitivity of touch	★★★★☆ (2.25)		★★★★☆ (2.25)		★★★★☆ (2)		★★★★☆ (2)	
	Colour of glove does not impede clinical activity	★★★★☆ (2)		★★★★☆ (2)		★★★★☆ (2)		★★★★☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	ROBINSON HEALTHCARE LIMITED	Supplier	MEDICARE PRODUCTS LTD	Supplier	HEALTHLINE	Supplier	HEALTHLINE
NITRILE NON-STERILE EXAMINATION GLOVES									
		Brand	Nytraguard Chemo Pure	Brand	Nitrex Accelerator Free	Brand	Premier Prestige	Brand	Healthline
		MPC	9542	MPC	GN06m	MPC	P2761AF6	MPC	GNE100/M
		NPC	FTE1844	NPC	FTE1084	NPC	FTE1286	NPC	FTG183
		Newton	6N	Newton	6N	Newton	6N	Newton	6N
		EN374 tested	Yes	EN374 tested	No information provided	EN374 tested	Yes	EN374 tested	Yes
		Unit of Issue	100	Unit of Issue	200	Unit of Issue	150	Unit of Issue	100
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✗		✗		✗	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★★☆ (2.25)	
	User able to maintain dexterity and grip when wet	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★★☆ (2.25)	
	User retains sensitivity of touch	★★★★☆ (2.25)		★★★☆☆ (2)		★★★☆☆ (2)		★★★★☆ (2.25)	
	Colour of glove does not impede clinical activity	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	




		Supplier	HEALTHLINE	Supplier	ROBINSON HEALTHCARE LIMITED	Supplier	UNIGLOVES (UK) LIMITED	Supplier	SHERMOND(DIV BUNZL RET SUP)
NON-STERILE VINYL EXAMINATION GLOVES									
		Brand	Handsafeglove	Brand	Vinoguard Premium	Brand	Unicare Premium Vinyl	Brand	Shermond Premier
		MPC	Gn65/m	MPC	9587	MPC	GU0063 (M)	MPC	P5855
		NPC	FTG286	NPC	FTE1867	NPC	FTE2348	NPC	FTG051
		Unit of Issue	100	Unit of Issue	1000	Unit of Issue	1000	Unit of Issue	100
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✓		✗		✗	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★☆☆ (1.75)		★★★☆☆ (1.75)		★★★☆☆ (2)		★★★☆☆ (2)	
	User able to maintain dexterity and grip when wet	★★★☆☆ (1.5)		★★★☆☆ (1.75)		★★★☆☆ (2)		★★★☆☆ (2)	
	User retains sensitivity of touch	★★★☆☆ (2)		★★★☆☆ (1.5)		★★★☆☆ (1.75)		★★★☆☆ (1.75)	
	Colour of glove does not impede clinical activity	★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)		★★★☆☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	





		Supplier	HEALTHLINE	Supplier	SHERMOND(DIV BUNZL RET SUP)	Supplier	ROBINSON HEALTHCARE LIMITED	Supplier	SUPERMAX HEALTHCARE LIMITED
NON-STERILE VINYL EXAMINATION GLOVES								Awaiting photograph	
		Brand	Healthline	Brand	Premier Prestige	Brand	Vinoguard Synthetic	Brand	Aurelia Delight PF - M
		MPC	Gn66/m	MPC	P7311	MPC	9592	MPC	38227
		NPC	FTG228	NPC	FTG189	NPC	FTE1872	NPC	FTE2085
		Unit of Issue	100	Unit of Issue	100	Unit of Issue	1000	Unit of Issue	1000
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✗		✗		✗	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2)		★★★★ (2)		★★★★ (2.25)		★★★★ (2)	
	User able to maintain dexterity and grip when wet	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (1.5)	
	User retains sensitivity of touch	★★★★ (1.75)		★★★★ (2)		★★★★ (2)		★★★★ (2)	
	Colour of glove does not impede clinical activity	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	





		Supplier	SEMPERIT	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	SHERMOND(DIV BUNZL RET SUP)	Supplier	FANNIN UK LTD
NON-STERILE VINYL EXAMINATION GLOVES									
		Brand	Sempermed	Brand	Aurelia Unique - Stretch Vinyl - M	Brand	Healthgard	Brand	Ice Clear Vinyl
		MPC	825781035	MPC	39227	MPC	PSCR0106	MPC	D1302 -11
		NPC	FTG197	NPC	FTE2089	NPC	FTG354	NPC	FTG443
		Unit of Issue	100	Unit of Issue	1000	Unit of Issue	100	Unit of Issue	100
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✗		✗		✓	
	Able to don easily on dry hands	✓		✓		✗		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2)		★★★★ (2)		★★★★ (1.5)		★★★★ (2)	
	User able to maintain dexterity and grip when wet	★★★★ (2)		★★★★ (2)		★★★★ (1.5)		★★★★ (1.75)	
	User retains sensitivity of touch	★★★★ (2)		★★★★ (1.75)		★★★★ (1.5)		★★★★ (2)	
	Colour of glove does not impede clinical activity	★★★★ (2.25)		★★★★ (2)		★★★★ (2)		★★★★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	





		Supplier	ROBINSON HEALTHCARE LIMITED	Supplier	HEALTHLINE	Supplier	FANNIN UK LTD	Supplier	SHERMOND(DIV BUNZL RET SUP)
NON-STERILE VINYL EXAMINATION GLOVES									
		Brand	Vinoguard Essentials	Brand	Healthline	Brand	Polar stretch soft white pigmented vinyl	Brand	Healthgard
		MPC	9582	MPC	GN130/02	MPC	D1202 -11	MPC	PSCR0162
		NPC	FTE1863	NPC	FTG147	NPC	FTE1527	NPC	FTG471
		Unit of Issue	1000	Unit of Issue	100	Unit of Issue	100	Unit of Issue	100
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✓		✗		✗		✗	
	Able to don easily on dry hands	✗		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★ ★ ★ ☆ (2)		★ ★ ★ ☆ (1.75)		★ ★ ★ ☆ (1.5)		★ ★ ★ ☆ (2)	
	User able to maintain dexterity and grip when wet	★ ★ ★ ☆ (1.75)		★ ★ ★ ☆ (1.5)		★ ★ ★ ☆ (1.5)		★ ★ ★ ☆ (1.75)	
	User retains sensitivity of touch	★ ★ ★ ☆ (1.75)		★ ★ ★ ☆ (1.75)		★ ★ ★ ☆ (1.5)		★ ★ ★ ☆ (2)	
	Colour of glove does not impede clinical activity	★ ★ ★ ☆ (2)		★ ★ ★ ☆ (1.75)		★ ★ ★ ☆ (2)		★ ★ ★ ☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	HEALTHLINE	Supplier	MEDICARE PRODUCTS LTD	Supplier	HEALTHLINE	Supplier	HEALTHLINE
NON-STERILE VINYL EXAMINATION GLOVES				Awaiting photograph				Awaiting photograph	
		Brand	Healthline	Brand	Medisafe	Brand	Healthline	Brand	Healthline
		MPC	GN661/M	MPC	GV05M	MPC	GN100/02	MPC	GN621
		NPC	FTG233	NPC	FTE1306	NPC	FTG128	NPC	FTG029
		Unit of Issue	100	Unit of Issue	100	Unit of Issue	100	Unit of Issue	100
		Score		Score		Score		Score	
Packaging and Storage	Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English		✓		✓		✓	
Opening and preparation for clinical use	Donning and fit	Opening large enough to remove gloves with ease		✓		✓		✓	
		Opening strong enough to resist tearing from repeated removal		✓		✓		✓	
		Product should fit in a dispenser		✓		✓		✓	
		Able to remove one glove at a time without wastage (eg multiple gloves falling out)		✗		✗		✗	
		Able to don easily on dry hands		✓		✓		✗	
Clinical Use	Clinical use, colour and material	User able to maintain dexterity and grip when dry		★★★★☆ (1.5)		★★★★☆ (2)		★★★★☆ (1.75)	
		User able to maintain dexterity and grip when wet		★★★★☆ (1.5)		★★★★☆ (2)		★★★★☆ (1.25)	
		User retains sensitivity of touch		★★★★☆ (1.25)		★★★★☆ (2)		★★★★☆ (1.25)	
		Colour of glove does not impede clinical activity		★★★★☆ (1.75)		★★★★☆ (2)		★★★★☆ (2)	
Safe disposal in the clinical environment	Removal and disposal	Easy to remove		✓		✓		✓	
		Can be disposed of by standard routes		✓		✓		✓	

		Supplier	MEDICARE PRODUCTS LTD	Supplier	SEMPERIT	Supplier	Hpc Healthline Uk Ltd
NON-STERILE VINYL EXAMINATION GLOVES							
		Brand	Syntex	Brand	Semperit Syn-Stretch	Brand	Healthline
		MPC	GS05M	MPC	825793235	MPC	GN63M
		NPC	FTE1472	NPC	FTG237	NPC	FTG554
		Unit of Issue	100	Unit of Issue	100	Unit of Issue	100
		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✓		✗		✓	
	Able to don easily on dry hands	✗		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2)		★★★★ (2)		★★★★ (1.4)	
	User able to maintain dexterity and grip when wet	★★★★ (2)		★★★★ (2)		★★★★ (1.8)	
	User retains sensitivity of touch	★★★★ (1.5)		★★★★ (2)		★★★★ (1.4)	
	Colour of glove does not impede clinical activity	★★★★ (1.75)		★★★★ (2)		★★★★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓	

		Supplier	HEALTHLINE	Supplier	SHERMOND(DIV BUNZL RET SUP)	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	UNIGLOVES (UK) LIMITED
NON-STERILE LATEX EXAMINATION GLOVES									
		Brand	Healthline	Brand	Shermond Premier Prestige	Brand	Aurelia Vibrant - M	Brand	Unicare Premium Latex
		MPC	Gn33/m	MPC	5355	MPC	98227	MPC	GU0013 (M)
		NPC	FTG118	NPC	FTE580	NPC	FTE2022	NPC	FTE2288
		Unit of Issue	100	Unit of Issue	25	Unit of Issue	1000	Unit of Issue	1000
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✗		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✓		✗		✓		✓	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★☆ (2.5)		★★★★☆ (1.75)		★★★★☆ (2.25)		★★★★☆ (2.25)	
	User able to maintain dexterity and grip when wet	★★★★☆ (2.5)		★★★★☆ (2)		★★★★☆ (2)		★★★★☆ (2.25)	
	User retains sensitivity of touch	★★★★☆ (2.25)		★★★★☆ (2)		★★★★☆ (2.25)		★★★★☆ (2.5)	
	Colour of glove does not impede clinical activity	★★★★☆ (2)		★★★★☆ (1.75)		★★★★☆ (2)		★★★★☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	SEMPERIT	Supplier	SHERMOND(DIV BUNZL RET SUP)
NON-STERILE LATEX EXAMINATION GLOVES									
		Brand	Aurelia Distinct - M	Brand	Aurelia Refresh - M	Brand	Sempermed	Brand	Premier Prestige
		MPC	29227	MPC	99227	MPC	823781735	MPC	P6602
		NPC	FTE2012	NPC	FTE2027	NPC	FTG090	NPC	FTE721
		Unit of Issue	1000	Unit of Issue	1000	Unit of Issue	100	Unit of Issue	100
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✗		✓		✓	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (3)		★★★★ (2.5)		★★★★ (2.25)		★★★★ (2.25)	
	User able to maintain dexterity and grip when wet	★★★★ (2.5)		★★★★ (2.5)		★★★★ (2.25)		★★★★ (2.25)	
	User retains sensitivity of touch	★★★★ (2.75)		★★★★ (2.5)		★★★★ (2.25)		★★★★ (2.25)	
	Colour of glove does not impede clinical activity	★★★★ (2.25)		★★★★ (2.25)		★★★★ (2)		★★★★ (2.25)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	HEALTHLINE	Supplier	HEALTHLINE	Supplier	MEDICARE PRODUCTS LTD	Supplier	FANNIN UK LTD
NON-STERILE LATEX EXAMINATION GLOVES									
		Brand	Healthline	Brand	Healthline	Brand	Bioplus	Brand	Dermagrip Ultra Latex
		MPC	GN30/m	MPC	Gn31/m	MPC	GL05M	MPC	D1002-20
		NPC	FTG086	NPC	FTE483	NPC		NPC	FTE1314
		Unit of Issue	100	Unit of Issue	100	Unit of Issue	100	Unit of Issue	100
		Score		Score		Score		Score	
Packaging and Storage	Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English		✓		✓		✓	
Opening and preparation for clinical use	Donning and fit	Opening large enough to remove gloves with ease		✓		✓		✓	
		Opening strong enough to resist tearing from repeated removal		✓		✓		✓	
		Product should fit in a dispenser		✓		✓		✓	
		Able to remove one glove at a time without wastage (eg multiple gloves falling out)		✗		✗		✓	
		Able to don easily on dry hands		✓		✓		✓	
Clinical Use	Clinical use, colour and material	User able to maintain dexterity and grip when dry		★★★★ (2.25)		★★★★ (2.25)		★★★★ (2.25)	
		User able to maintain dexterity and grip when wet		★★★★ (2.25)		★★★★ (2.25)		★★★★ (2.5)	
		User retains sensitivity of touch		★★★★ (2.25)		★★★★ (2.25)		★★★★ (2)	
		Colour of glove does not impede clinical activity		★★★★ (2.25)		★★★★ (2.25)		★★★★ (2.25)	
Safe disposal in the clinical environment	Removal and disposal	Easy to remove		✓		✓		✓	
		Can be disposed of by standard routes		✓		✓		✓	

		Supplier	SHERMOND(DIV BUNZL RET SUP)	Supplier	ROBINSON HEALTHCARE LIMITED	Supplier	SUPERMAX HEALTHCARE LIMITED	Supplier	UNIGLOVES (UK) LIMITED
NON-STERILE LATEX EXAMINATION GLOVES									
		Brand	Premier Protector	Brand	Cyraguard Essentials	Brand	Aurelia Velocity Original - M	Brand	Supergrip
		MPC	P7902	MPC	9502	MPC	88267	MPC	GM0023 (M)
		NPC	FTG365	NPC	FTE1855	NPC	FTE2017	NPC	FTE2293
		Unit of Issue	100	Unit of Issue	1000	Unit of Issue	1000	Unit of Issue	1000
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓		✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓		✓		✓	
	Product should fit in a dispenser	✓		✓		✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✓		✓		✓	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2.25)		★★★★ (2)		★★★★ (2.25)		★★★★ (2.25)	
	User able to maintain dexterity and grip when wet	★★★★ (2.25)		★★★★ (2)		★★★★ (2.25)		★★★★ (2)	
	User retains sensitivity of touch	★★★★ (2.5)		★★★★ (2.5)		★★★★ (2.5)		★★★★ (2.5)	
	Colour of glove does not impede clinical activity	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	FANNIN UK LTD	Supplier	ROBINSON HEALTHCARE LIMITED
NON-STERILE LATEX EXAMINATION GLOVES					
		Brand	Dermagrip	Brand	Cyraguard Premium
		MPC	D1402-14	MPC	9507
		NPC	FTE1796	NPC	FTE1859
		Unit of Issue	500	Unit of Issue	1000
		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓	
Opening and preparation for clinical use Donning and fit	Opening large enough to remove gloves with ease	✓		✓	
	Opening strong enough to resist tearing from repeated removal	✓		✓	
	Product should fit in a dispenser	✓		✓	
	Able to remove one glove at a time without wastage (eg multiple gloves falling out)	✗		✗	
	Able to don easily on dry hands	✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★☆ (2)		★★★★☆ (2.25)	
	User able to maintain dexterity and grip when wet	★★★★☆ (1.75)		★★★★☆ (2.25)	
	User retains sensitivity of touch	★★★★☆ (1.75)		★★★★☆ (2.5)	
	Colour of glove does not impede clinical activity	★★★★☆ (1.75)		★★★★☆ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓	
	Can be disposed of by standard routes	✓		✓	

		Supplier	SHERMOND(DIV BUNZL RET SUP)	Supplier	MEDICARE PRODUCTS LTD	Supplier	MEDICARE PRODUCTS LTD	Supplier	FANNIN UK LTD
NITRILE STERILE EXAMINATION GLOVES									
		Brand	Premier	Brand	Nitrex Sterile	Brand	Intouch V Synthetic	Brand	Dermagrip Sterile Nitrile
		MPC	P2761/S	MPC	GN35M	MPC	W02EM	MPC	D1502-82
		NPC	FTG498	NPC	FTE1531	NPC	FTG447	NPC	FTE1310
		Newton	6N	Newton	No info available	Newton	6N	Newton	6N
		EN374 tested	Y	EN374 tested	No info available	EN374 tested	No info available	EN374 tested	No info available
		Unit of Issue	200	Unit of Issue	50	Unit of Issue	50	Unit of Issue	50
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✗		✗		✓	
PACKAGING AND STORAGE Internal packet	The product type is clearly visible - pair of sterile exam gloves	✓		✓		✓		✓	
	Instructions for opening are clear	✓		✗		✗		✓	
	The glove material is clearly visible	✓		✓		✓		✓	
	The size is clearly visible	✓		✓		✓		✓	
	The expiry date, EN standards and CE marking are clearly visible	✓		✓		✓		✓	
	The information is in English and clearly explains what is in the packet	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Able to open packet without contaminating	✓		✓		✓		✓	
	Able to access the inner packet without contaminating	✓		✓		✓		✓	
	Able to open inner packet without contaminating	✓		✓		✓		✓	
	Able to identify L and R hand	✓		✓		✓		✓	
	Able to don easily on dry hands	✓		✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (1.75)		★★★★ (2)		★★★★ (2)		★★★★ (1.5)	
	User able to maintain dexterity and grip when wet	★★★★ (1.75)		★★★★ (2)		★★★★ (2)		★★★★ (1.5)	
	User retains sensitivity of touch	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (1.5)	
	Colour of glove does not impede clinical activity	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	SHERMOND(DIV BUNZL RET SUP)	Supplier	Handsafe	Supplier	Handsafe	Supplier	SEMPERIT
NITRILE STERILE EXAMINATION GLOVES								Awaiting Photograph	
		Brand	Premier Prestige	Brand	Handsafe	Brand	Handsafe	Brand	Sempercure
		MPC	P2761AF/S	MPC	GS6906	MPC	Sne100/02	MPC	826753727
		NPC	FTE798	NPC	FTG203	NPC	FTG250	NPC	FTG289
		Newton	9N	Newton	6N	Newton	No info available	Newton	9N
		EN374 tested	No info available	EN374 tested	Y	EN374 tested	No info available	EN374 tested	No info available
		Unit of Issue	50	Unit of Issue	50	Unit of Issue	50	Unit of Issue	50
		Score		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓	
PACKAGING AND STORAGE Internal packet	The product type is clearly visible - pair of sterile exam gloves	✓		✓		✓		✓	
	Instructions for opening are clear	✓		✓		✓		✓	
	The glove material is clearly visible	✓		✓		✓		✓	
	The size is clearly visible	✓		✓		✗		✓	
	The expiry date, EN standards & CE marking are clearly visible	✓		✗		✓		✓	
	The information is in English and clearly explains what is in the packet	✓		✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Able to open packet without contaminating	✓		✓		✓		✓	
	Able to access the inner packet without contaminating	✓		✓		✓		✓	
	Able to open inner packet without contaminating	✓		✓		✓		✓	
	Able to identify L and R hand	✓		✓		✓		✓	
	Able to don easily on dry hands	✓		✓		✗		✗	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2)		★★★★ (2)		★★★★ (1)		★★★★ (1.5)	
	User able to maintain dexterity and grip when wet	★★★★ (2)		★★★★ (2)		★★★★ (1.5)		★★★★ (1.5)	
	User retains sensitivity of touch	★★★★ (2)		★★★★ (2)		★★★★ (1.5)		★★★★ (1.5)	
	Colour of glove does not impede clinical activity	★★★★ (2)		★★★★ (2)		★★★★ (1)		★★★★ (1.5)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓		✓	

		Supplier	SHERMOND(DIV BUNZL RET SUP)
NITRILE STERILE ACCELERATOR FREE EXAMINATION GLOVES			
		Brand	Premier Prestige
		MPC	P2761AF/S
		NPC	FTE798
		Newton	9N
		EN374 tested	No information provided
		Unit of Issue	50
		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓	
PACKAGING AND STORAGE Internal packet	The product type is clearly visible - pair of sterile exam gloves	✓	
	Instructions for opening are clear	✓	
	The glove material is clearly visible	✓	
	The size is clearly visible	✓	
	The expiry date, EN standards and CE marking are clearly visible	✓	
	The information is in English and clearly explains what is in the packet	✓	
Opening and preparation for clinical use Donning and fit	Able to open packet without contaminating	✓	
	Able to access the inner packet without contaminating	✓	
	Able to open inner packet without contaminating	✓	
	Able to identify L and R hand	✓	
	Able to don easily on dry hands	✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★ ★ ★ (2)	
	User able to maintain dexterity and grip when wet	★ ★ ★ (2)	
	User retains sensitivity of touch	★ ★ ★ (2)	
	Colour of glove does not impede clinical activity	★ ★ ★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓	
	Can be disposed of by standard routes	✓	

		Supplier	HANDSAFE	Supplier	SHERMOND (DIV BUNZL RET SUP)
VINYL STERILE EXAMINATION GLOVES					
		Brand	Handsafe	Brand	Shermond Premier SV Sterile
		MPC	GS124/02	MPC	P4302
		NPC	FTE810	NPC	FTG213
		Unit of Issue	50	Unit of Issue	50
		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓	
PACKAGING AND STORAGE Internal packet	The product type is clearly visible - pair of sterile exam gloves	✓		✓	
	Instructions for opening are clear	✓		✓	
	The glove material is clearly visible	✓		✓	
	The size is clearly visible	✓		✓	
	The expiry date, EN standards and CE marking are clearly visible	✓		✓	
	The information is in English and clearly explains what is in the packet	✓		✓	
Opening and preparation for clinical use Donning and fit	Able to open packet without contaminating	✓		✓	
	Able to access the inner packet without contaminating	✓		✓	
	Able to open inner packet without contaminating	✓		✓	
	Able to identify L and R hand	✓		✓	
	Able to don easily on dry hands	✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★ ★ ★ (1.75)		★ ★ ★ (2)	
	User able to maintain dexterity and grip when wet	★ ★ ★ (1.75)		★ ★ ★ (1.75)	
	User retains sensitivity of touch	★ ★ ★ (1.75)		★ ★ ★ (1.75)	
	Colour of glove does not impede clinical activity	★ ★ ★ (2)		★ ★ ★ (2)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓	
	Can be disposed of by standard routes	✓		✓	

		Supplier	SEMPERIT	Supplier	Handsafe	Supplier	MEDICARE PRODUCTS LTD	Supplier	SHERMOND(DIV BUNZL RET SUP)
LATEX STERILE EXAMINATION GLOVES									
		Brand	Sempermed	Brand	Handsafe	Brand	Intouch V Natural	Brand	Shermond Premier
		MPC	RT826759727	MPC	Gs21/m	MPC	INTW80EM	MPC	P3302
		NPC	FTG123	NPC	FTG023	NPC	FTG482	NPC	FTE883
		Unit of Issue	40	Unit of Issue	50	Unit of Issue	300	Unit of Issue	50
		Score		Score		Score		Score	
Packaging and Storage	Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓		✓
PACKAGING AND STORAGE	Internal packet	The product type is clearly visible - pair of sterile exam gloves	✓		✓		✓		✓
		Instructions for opening are clear	✗		✓		✗		✓
		The glove material is clearly visible	✓		✓		✓		✓
		The size is clearly visible	✓		✓		✓		✓
		The expiry date, EN standards and CE marking are clearly visible	✗		✓		✗		✓
		The information is in English and clearly explains what is in the packet	✓		✓		✓		✓
Opening and preparation for clinical use	Donning and fit	Able to open packet without contaminating	✓		✓		✓		✓
		Able to access the inner packet without contaminating	✓		✓		✓		✓
		Able to open inner packet without contaminating	✓		✓		✓		✓
		Able to identify L and R hand	✓		✓		✓		✓
		Able to don easily on dry hands	✓		✓		✓		✓
Clinical Use	Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★★ (2.25)		★★★★ (2)		★★★★ (2)		★★★★ (2.25)
		User able to maintain dexterity and grip when wet	★★★★ (2.25)		★★★★ (2)		★★★★ (2)		★★★★ (2.25)
		User retains sensitivity of touch	★★★★ (2.25)		★★★★ (2.25)		★★★★ (2.25)		★★★★ (2.25)
		Colour of glove does not impede clinical activity	★★★★ (2)		★★★★ (2)		★★★★ (2)		★★★★ (2)
Safe disposal in the clinical environment	Removal and disposal	Easy to remove	✓		✓		✓		✓
		Can be disposed of by standard routes	✓		✓		✓		✓

		Supplier	Handsafe	Supplier	Fannin	Supplier	MEDICARE PRODUCTS LTD
LATEX STERILE EXAMINATION GLOVES							
		Brand	Handsafe	Brand	Dermagrip	Brand	Neotex
		MPC	GS680	MPC	D1502-13a	MPC	GS36M
		NPC	FTG014	NPC	FTE1318	NPC	FTG417
		Unit of Issue	50	Unit of Issue	50	Unit of Issue	200
		Score		Score		Score	
Packaging and Storage Box	The product type, size, expiry date, material, EN standards, lot numbers and CE markings are clearly visible and in English	✓		✓		✓	
PACKAGING AND STORAGE Internal packet	The product type is clearly visible - pair of sterile exam gloves	✓		✓		✓	
	Instructions for opening are clear	✓		✓		✗	
	The glove material is clearly visible	✓		✓		✗	
	The size is clearly visible	✓		✓		✓	
	The expiry date, EN standards and CE marking are clearly visible	✗		✓		✗	
	The information is in English and clearly explains what is in the packet	✓		✓		✓	
Opening and preparation for clinical use Donning and fit	Able to open packet without contaminating	✓		✓		✗	
	Able to access the inner packet without contaminating	✓		✓		✗	
	Able to open inner packet without contaminating	✓		✓		✗	
	Able to identify L and R hand	✓		✓		✓	
	Able to don easily on dry hands	✓		✓		✓	
Clinical Use Clinical use, colour and material	User able to maintain dexterity and grip when dry	★★★☆☆ (2.00)		★★★☆☆ (2.00)		★★★☆☆ (1.75)	
	User able to maintain dexterity and grip when wet	★★★★☆ (2.25)		★★★☆☆ (2.00)		★★★★☆ (2.00)	
	User retains sensitivity of touch	★★★★☆ (2.25)		★★★★☆ (2.00)		★★★☆☆ (1.50)	
	Colour of glove does not impede clinical activity	★★★☆☆ (2.00)		★★★☆☆ (2.00)		★★★☆☆ (1.75)	
Safe disposal in the clinical environment Removal and disposal	Easy to remove	✓		✓		✓	
	Can be disposed of by standard routes	✓		✓		✓	