










PROCEDURE UNDERPADS		Supplier	Abena UK Ltd			Abena UK Ltd			Attends Ltd					Attends Ltd	
															
NHS CET Scoring Matrix		Brand	Abri Soft Eco			Abri Soft Basic			Attends Cover Dri - Plus					Attends Cover Dri - Super	
NHS CET Product Assessment Cycle		MPC	254116	254117	254118	4116	4117	4118	203910	203934	203972	204016	204030	203958	203996
Assessment Criteria		NPC							VMU331	VMU022	VMU026		VMU031	VMU032	VMU033
Description		Unit of Issue	Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm & 60x90cm			Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm & 60x90cm			Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm, 60x90cm, 80x90cm & 170x90cm					Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x60 & 60x90cm	
			240	240	120	240	240	120	50	50	50	30	50	50	50
Packaging	The packaging has an easily identified description of correct use and size. All necessary detail such as Lot Number and Manufacturer are easy to find	Scores	★★★★ (2.80)			★★★★ (2.80)			★★★ (1.60)					★★★ (1.60)	
	The packaging can be easily opened.	Scores	★★★★ (2.60)			★★★ (1.60)			★★★★ (2.80)					★★★★ (2.40)	
Opening and Preparation	The pads are dispensed individually	Scores	★★★★ (3.00)			★★★★ (2.60)			★★★★ (2.80)					★★★★ (2.80)	
	The packaging remains robust to support the products in a useable condition until empty	Scores	★★★★ (2.60)			★★★ (2.20)			★★★★ (3.00)					★★★★ (3.00)	
Clinical Use Examples	The pad in its dry state is comfortable for a procedure	Scores	★★★★ (2.60)			★★★★ (2.40)			★★★★ (2.60)					★★★★ (2.60)	
	The colour of the absorbent surface does not impede clinical decision making	Scores	★★★★ (2.80)			★★★★ (2.80)			★★★★ (2.80)					★★★★ (2.80)	
Light absorption required Examinations, insertion of lines, skin prep, dressing & wound procedures, catheter procedures, oral secretions, bottom end procedures, hygiene needs, formed faecal incontinence and endoscopy		Sizes	60x40cm	60x60cm	60x90cm	60x40cm	60x60cm	60x90cm	60x40cm	60x60cm	60x90cm	80x90cm	170x90cm	60x60cm	60x90cm
	Fluid management - light absorbency	Scores	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Moderate absorption required Maternity-postpartum, theatre blood loss, wound wash out and oedematous patients / limbs	Fluid management - moderate absorbency	Scores	✗	✓	✓	✗	✓	✓	✗	✓	✓	✓	✓	✓	✓
	Fluid management - high absorbency	Scores	✗	✗	✗	✗	✗	✓	✗	✗	✗	✓	✓	✗	✓
High absorption required Loose faecal incontinence, Maternity-intrapartum, theatre wash-out, theatre trauma blood loss and loss of skin barrier e.g. burns fluid seep	Surface of the pad to remain dry when pad has been moistened	Scores	★★ (0.00)**			★★ (0.00)**			★★ (0.00)**					★★ (0.20)**	
	The pad to remain waterproof to protect the surface underneath it	Scores	★★★★ (2.00)**			★★★★ (2.00)**			★★★★ (2.00)**					★★★★ (2.00)**	
Disposal	Retention of liquid within the pad	Scores	★★★ (1.20)			★★★★ (2.00)			★★★★ (2.00)					★★★★ (2.00)	
	Pad can be disposed safely	Scores	★★★ (1.40)**			★★★★ (2.00)**			★★★★ (2.00)**					★★★★ (2.00)**	


**Maximum number of stars attainable 2

PROCEDURE UNDERPADS		Supplier	Drylock Technologies		Drylock Technologies		Drylock Technologies		Drylock Technologies						
															
NHS CET Scoring Matrix		Brand	Dailee combed premium fix		Dailee combed Normal		Dailee combed premium air		Dailee combed Plus						
Score		Meaning		MPC	626205		622202	624204	626204	626402		622301	624302	626303	
0		This does not meet the criteria		NPC	VMU403		VMU397	VMU399	VMU401	VMU404		VMU398	VMU400	VMU402	
1		The partially meets the criteria		Description		Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Size 60x90cm		Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm, 60x90cm		Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Size 60x90cm		Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Size 60x40cm, 60x60cm & 60x90cm			
2		This meets the criteria		Unit of Issue		25		25	25	25	30		25	25	25
3		This exceeds the criteria		Scores		Pre-production sample only available at time of evaluation		Pre-production sample only available at time of evaluation		Pre-production sample only available at time of evaluation		Pre-production sample only available at time of evaluation			
Packaging	The packaging has an easily identified description of correct use and size. All necessary detail such as Lot Number and Manufacturer are easy to find	Scores	Pre-production sample only available at time of evaluation		Pre-production sample only available at time of evaluation		Pre-production sample only available at time of evaluation		Pre-production sample only available at time of evaluation		Pre-production sample only available at time of evaluation				
	The packaging can be easily opened.	Scores	★ ★ ★ (1.40)		★ ★ ★ (1.40)		★ ★ ★ (1.40)		★ ★ ★ (1.40)		★ ★ ★ (1.40)				
Opening and Preparation	The pads are dispensed individually	Scores	★ ★ ★ (1.60)		★ ★ ★ (1.80)		★ ★ ★ (1.80)		★ ★ ★ (1.80)		★ ★ ★ (1.60)				
	The packaging remains robust to support the products in a useable condition until empty	Scores	★ ★ ★ (1.40)		★ ★ ★ (1.60)		★ ★ ★ (1.60)		★ ★ ★ (1.60)		★ ★ ★ (1.40)				
Clinical Use	The pad in its dry state is comfortable for a procedure	Scores	★ ★ ★ ★ (3.00)		★ ★ ★ ★ (3.00)		★ ★ ★ ★ (2.80)		★ ★ ★ ★ (3.00)		★ ★ ★ ★ (3.00)				
	The colour of the absorbent surface does not impede clinical decision making	Scores	★ ★ ★ ★ (3.00)		★ ★ ★ ★ (2.80)		★ ★ ★ ★ (2.80)		★ ★ ★ ★ (3.00)		★ ★ ★ ★ (3.00)				
Examples	Light absorption required	Sizes	60x90cm		60x40cm	60x60cm	60x90cm	60x90cm		60x40cm	60x60cm	60x90cm			
	Examinations, insertion of lines, skin prep, dressing & wound procedures, catheter procedures, oral secretions, bottom end procedures, hygiene needs, formed faecal incontinence and endoscopy	Fluid management - light absorbency	Scores	✓		✓	✓	✓	✓		✓	✓	✓		
Moderate absorption required	Fluid management - moderate absorbency	Scores	✓		✗	✓	✓	✓		✗	✓	✓			
	Maternity-postpartum, theatre blood loss, wound wash out and oedematous patients / limbs	Fluid management - high absorbency	Scores	✗		✗	✗	✗	✓		✗	✗	✓		
High absorption required	Surface of the pad to remain dry when pad has been moistened	Scores	★ ★ (0.00) **		★ ★ (0.20) **		★ ★ (0.00) **		★ ★ (0.40) **		★ ★ (0.40) **				
	Loose faecal incontinence, Maternity-intrapartum, theatre wash-out, theatre trauma blood loss and loss of skin barrier e.g. burns fluid seep	The pad to remain waterproof to protect the surface underneath it	Scores	★ ★ ★ (2.00) **		★ ★ ★ (2.00) **		★ ★ ★ (2.00) **		★ ★ ★ (2.00) **		★ ★ ★ (2.00) **			
Disposal	Retention of liquid within the pad	Scores	★ ★ ★ (2.00)		★ ★ ★ (2.00)		★ ★ ★ (1.80)		★ ★ ★ (2.00)		★ ★ ★ (2.00)				
	Pad can be disposed safely	Scores	★ ★ ★ (1.80) **		★ ★ ★ (2.00) **		★ ★ ★ (1.80) **		★ ★ ★ (2.00) **		★ ★ ★ (2.00) **				




**Maximum number of stars attainable 2

PROCEDURE UNDERPADS		Supplier	Fourstones Paper Mill Co Ltd				
							
NHS CET Scoring Matrix		Brand	Warden BFM				
		MPC	6415DB	6479DB	VMU096	6090DB	
		NPC	VMU260	VMU261	VMU096	VMU262	
NHS CET Product Assessment Cycle		Assessment Criteria	Description	Disposable procedure underpad with Cellulose wadding Recycled Pulp 4 way sealed & waterproof backing. Size 40x60cm, 60x55cm, 60x75cm & 60x90cm			
		Unit of Issue	100	100	100	100	
Packaging	The packaging has an easily identified description of correct use and size. All necessary detail such as Lot Number and Manufacturer are easy to find	Scores	★☆☆ (1.20)				
	The packaging can be easily opened.	Scores	★★★☆☆ (1.80)				
Opening and Preparation	The pads are dispensed individually	Scores	★★★☆☆ (2.20)				
	The packaging remains robust to support the products in a useable condition until empty	Scores	★★★★☆ (2.40)				
Clinical Use Examples Light absorption required Examinations, insertion of lines, skin prep, dressing & wound procedures, catheter procedures, oral secretions, bottom end procedures, hygiene needs, formed faecal incontinence and endoscopy Moderate absorption required Maternity-postpartum, theatre blood loss, wound wash out and oedematous patients / limbs High absorption required Loose faecal incontinence, Maternity-intrapartum, theatre wash-out, theatre trauma blood loss and loss of skin barrier e.g. burns fluid seep	The pad in its dry state is comfortable for a procedure	Scores	★★★☆☆ (1.40)				
	The colour of the absorbent surface does not impede clinical decision making	Scores	☆☆☆☆ (0.20)				
		Sizes	60x40cm	60x60cm	60x75cm	60x90cm	
	Fluid management - light absorbency	Scores	✓	✓	✓	✓	
	Fluid management - moderate absorbency	Scores	✗	✗	✓	✓	
	Fluid management - high absorbency	Scores	✗	✗	✗	✗	
	Surface of the pad to remain dry when pad has been moistened	Scores	☆☆☆☆ (0.00) **				
	The pad to remain waterproof to protect the surface underneath it	Scores	★★★☆☆ (1.60) **				
	Retention of liquid within the pad	Scores	☆☆☆☆ (0.60)				
Disposal	Pad can be disposed safely	Scores	☆☆☆☆ (0.80) **				



**Maximum number of stars attainable 2

PROCEDURE UNDERPADS		Supplier	Mahr Impex UK Ltd			
						
Score	Meaning	Brand	ADA COMFORT			
0	This does not meet the criteria	MPC	79000150604090	79000150606090	79000200000000	79000150607590
1	The partially meets the criteria	NPC	VMU423	VMU424		VMU425
2	This meets the criteria	Description	Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm, 60x90cm & 170x90cm			
3	This exceeds the criteria	Unit of Issue	240	150	150	150
NHS CET Product Assessment Cycle	Assessment Criteria	Description	Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm, 60x90cm & 170x90cm			
Packaging	The packaging has an easily identified description of correct use and size. All necessary detail such as Lot Number and Manufacturer are easy to find	Scores	★ ★ ★ (1.60)			
	The packaging can be easily opened.	Scores	★ ★ ★ (1.80)			
Opening and Preparation	The pads are dispensed individually	Scores	★ ★ ★ (2.20)			
	The packaging remains robust to support the products in a useable condition until empty	Scores	★ ★ ★ (1.60)			
Clinical Use Examples Light absorption required Examinations, insertion of lines, skin prep, dressing & wound procedures, catheter procedures, oral secretions, bottom end procedures, hygiene needs, formed faecal incontinence and endoscopy Moderate absorption required Maternity-postpartum, theatre blood loss, wound wash out and oedematous patients / limbs High absorption required Loose faecal incontinence, Maternity-intrapartum, theatre wash-out, theatre trauma blood loss and loss of skin barrier e.g. burns fluid seep	The pad in its dry state is comfortable for a procedure	Scores	★ ★ ★ (1.80)			
	The colour of the absorbent surface does not impede clinical decision making	Scores	★ ★ ★ (2.80)			
		Sizes	60x40cm	60x60cm	60x75cm	60x90cm
	Fluid management - light absorbency	Scores	✓	✓	✓	✓
	Fluid management - moderate absorbency	Scores	✗	✓	✓	✓
	Fluid management - high absorbency	Scores	✗	✗	✗	✗
	Surface of the pad to remain dry when pad has been moistened	Scores	★ ★ (0.00) **			
	The pad to remain waterproof to protect the surface underneath it	Scores	★ ★ ★ (2.00) **			
	Retention of liquid within the pad	Scores	★ ★ ★ (1.80)			
Disposal	Pad can be disposed safely	Scores	★ ★ ★ (1.80) **			



**Maximum number of stars attainable 2

PROCEDURE UNDERPADS		Supplier	Medi-Inn Ltd			Medi-Inn Ltd				Medi-Inn Ltd			
													
NHS CET Scoring Matrix		Brand	Medi-Inn Plus			Medi-Inn Super				Medi-Inn Cellulose			
NHS CET Product Assessment Cycle		MPC	N12500	N12600	N12700	N11069	N11071	N11703	N11000	N12000	N12001	N12002	N12003
Assessment Criteria		NPC	VMU416	VMU418	VMU420	VMU417	VMU419	VMU421	VMU422	VMU412	VMU413	VMU414	VMU415
Description		Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm & 60x90cm			Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm, 60x90cm & 170x90cm				Disposable procedure underpad with Cellulose wadding Recycled Pulp 4 way sealed & waterproof backing. Size 40x60cm, 60x60cm, 60x75cm & 60x90cm				
Unit of Issue		300 150 150			300 150 150 100				200 150 150 100				
Packaging	The packaging has an easily identified description of correct use and size. All necessary detail such as Lot Number and Manufacturer are easy to find	Scores	Pre-production sample only available at time of evaluation			Pre-production sample only available at time of evaluation				Pre-production sample only available at time of evaluation			
	The packaging can be easily opened.	Scores	★ ★ ★ (1.80)			★ ★ ★ (2.20)				★ ★ ★ (1.80)			
Opening and Preparation	The pads are dispensed individually	Scores	★ ★ ★ (2.80)			★ ★ ★ (2.20)				★ ★ ★ (2.60)			
	The packaging remains robust to support the products in a useable condition until empty	Scores	★ ★ ★ (2.40)			★ ★ ★ (2.00)				★ ★ ★ (2.40)			
Clinical Use	The pad in its dry state is comfortable for a procedure	Scores	★ ★ ★ (1.80)			★ ★ ★ (1.80)				★ ★ ★ (1.80)			
	The colour of the absorbent surface does not impede clinical decision making	Scores	★ ★ ★ (2.20)			★ ★ ★ (3.00)				★ ★ ★ (2.40)			
Examples		Sizes	60x40cm	60x60cm	60x90cm	60x40cm	60x60cm	60x90cm	170x90cm	60x40cm	60x60cm	60x75cm	60x90cm
Light absorption required	Examinations, insertion of lines, skin prep, dressing & wound procedures, catheter procedures, oral secretions, bottom end procedures, hygiene needs, formed faecal incontinence and endoscopy	Scores	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
Moderate absorption required	Fluid management - moderate absorbency	Scores	✗	✓	✓	✗	✓	✓	✓	✗	✗	✓	✓
	Fluid management - high absorbency	Scores	✗	✗	✓	✗	✗	✓	✓	✗	✗	✗	✗
High absorption required	Maternity-postpartum, theatre blood loss, wound wash out and oedematous patients / limbs	Scores	★ ★ (0.00) **			★ ★ (0.40) **				★ ★ (0.00) **			
	The pad to remain waterproof to protect the surface underneath it	Scores	★ ★ (2.00) **			★ ★ (2.00) **				★ ★ (2.00) **			
	Retention of liquid within the pad	Scores	★ ★ ★ (2.00)			★ ★ ★ (2.00)				★ ★ ★ (1.60)			
Disposal	Pad can be disposed safely	Scores	★ ★ (1.80) **			★ ★ (2.00) **				★ ★ (1.40) **			




**Maximum number of stars attainable 2

PROCEDURE UNDERPADS		Supplier	Ontex Healthcare UK Ltd			Ontex Healthcare UK Ltd						
												
NHS CET Scoring Matrix		Brand	ID Protect Expert Plus			ID protect expert Super						
Score		Meaning		MPC	5800460300	5800660300	5800960300	5800475300	5800675300	5800775300	5800975300	5800075200
0		This does not meet the criteria		NPC	VMU366	CFP1741	VMU363	VMU388	VMU365	VMU364	VMU361	CFP1734
1		The partially meets the criteria		Description	Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm & 60x90cm			Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm, 60x75cm, 60x90cm, 90x180cm				
2		This meets the criteria		Unit of Issue	30	30	30	270	30	30	30	20
3		This exceeds the criteria										
NHS CET Product Assessment Cycle	Assessment Criteria	Scores										
Packaging	The packaging has an easily identified description of correct use and size. All necessary detail such as Lot Number and Manufacturer are easy to find	Scores	★★★★ (2.60)			★★★★ (2.60)						
	The packaging can be easily opened.	Scores	★★★★ (2.40)			★★★★ (2.60)						
Opening and Preparation	The pads are dispensed individually	Scores	★★★☆☆ (2.00)			★★★☆☆ (2.20)						
	The packaging remains robust to support the products in a useable condition until empty	Scores	★★★☆☆ (2.20)			★★★☆☆ (2.20)						
Clinical Use Examples Light absorption required Examinations, insertion of lines, skin prep, dressing & wound procedures, catheter procedures, oral secretions, bottom end procedures, hygiene needs, formed faecal incontinence and endoscopy Moderate absorption required Maternity-postpartum, theatre blood loss, wound wash out and oedematous patients / limbs High absorption required Loose faecal incontinence, Maternity-intrapartum, theatre wash-out, theatre trauma blood loss and loss of skin barrier e.g. burns fluid seep	The pad in its dry state is comfortable for a procedure	Scores	★★★★ (2.60)			★★★★ (2.40)						
	The colour of the absorbent surface does not impede clinical decision making	Scores	★★★★ (2.80)			★★★★ (2.80)						
		Sizes	60x40cm	60x60cm	60x90cm	60x40cm	60x60cm	60x75cm	60x90cm	90x180cm		
	Fluid management - light absorbency	Scores	✓	✓	✓	✓	✓	✓	✓	✓		
	Fluid management - moderate absorbency	Scores	✗	✓	✓	✗	✓	✓	✓	✓		
	Fluid management - high absorbency	Scores	✗	✗	✗	✗	✗	✗	✓	✓		
	Surface of the pad to remain dry when pad has been moistened	Scores	☆☆☆☆ (0.00) **			☆☆☆☆ (0.00) **						
	The pad to remain waterproof to protect the surface underneath it	Scores	★★★★ (2.00) **			★★★★ (2.00) **						
Retention of liquid within the pad	Scores	★★★☆☆ (1.60)			★★★☆☆ (2.00)							
Disposal	Pad can be disposed safely	Scores	★★★☆☆ (1.60) **			★★★☆☆ (2.00) **						

**Maximum number of stars attainable 2

PROCEDURE UNDERPADS		Supplier	Robinson Healthcare Limited				Robinson Healthcare Limited			
										
Score	Meaning	Brand	Readit				Readit			
0	This does not meet the criteria	MPC	5515	5516	5508	5425	5503	5504	5509	5426
1	The partially meets the criteria	NPC	VMU389	VMU391	VMU393	VMU395	VMU390	VMU392	VMU394	VMU396
2	This meets the criteria	Description	Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm 60x90cm & 170x90cm				Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm 60x90cm & 170x90cm			
3	This exceeds the criteria	Unit of Issue	100	100	100	100	100	100	100	100
NHS CET Scoring Matrix		NHS CET Product Assessment Cycle	Assessment Criteria							
Packaging	The packaging has an easily identified description of correct use and size. All necessary detail such as Lot Number and Manufacturer are easy to find	Scores	★★★★ (2.20)				★★★★ (2.20)			
	The packaging can be easily opened.	Scores	★★★★ (1.40)				★★★★ (1.40)			
Opening and Preparation	The pads are dispensed individually	Scores	★★★ (1.20)				★★★ (1.20)			
	The packaging remains robust to support the products in a useable condition until empty	Scores	★★★★ (1.60)				★★★★ (1.60)			
Clinical Use	The pad in its dry state is comfortable for a procedure	Scores	★★★★ (1.40)				★★★★ (1.40)			
	The colour of the absorbent surface does not impede clinical decision making	Scores	★★★★ (2.40)				★★★★ (2.60)			
Examples		Sizes	60x40cm	60x60cm	60x90cm	170x90cm	60x40cm	60x60cm	60x90cm	170x90cm
Light absorption required Examinations, insertion of lines, skin prep, dressing & wound procedures, catheter procedures, oral secretions, bottom end procedures, hygiene needs, formed faecal incontinence and endoscopy	Fluid management - light absorbency	Scores	✓	✓	✓	✓	✓	✓	✓	✓
	Fluid management - moderate absorbency	Scores	✗	✓	✓	✓	✗	✓	✓	✓
Moderate absorption required Maternity-postpartum, theatre blood loss, wound wash out and oedematous patients / limbs	Fluid management - high absorbency	Scores	✗	✗	✓	✓	✗	✗	✓	✓
	Surface of the pad to remain dry when pad has been moistened	Scores	★★ (0.00) **				★★ (0.20) **			
High absorption required Loose faecal incontinence, Maternity-intrapartum, theatre wash-out, theatre trauma blood loss and loss of skin barrier e.g. burns fluid seep	The pad to remain waterproof to protect the surface underneath it	Scores	★★★★ (1.80) **				★★★★ (2.00) **			
	Retention of liquid within the pad	Scores	★★★★ (1.80)				★★★★ (2.00)			
Disposal	Pad can be disposed safely	Scores	★★★ (1.60) **				★★★ (2.00) **			

**Maximum number of stars attainable 2

PROCEDURE UNDERPADS		Supplier	SCA UK Ltd		SCA UK Ltd					SCA UK Ltd	
											
NHS CET Scoring Matrix		Brand	TENA Bed Normal		TENA Bed Plus					TENA Bed Super	
NHS CET Product Assessment Cycle		MPC	770045	770047	770113	770100	770115	770104	771103	770202	770200
Assessment Criteria		NPC	VMU375	VMU376	VMU349	VMU309	VMU313	VMU310	VMU316	VMU315	VMU314
Description		Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x60cm & 60x90cm		Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x40cm, 60x60cm, 60x75cm, 60x90cm & 80x180cm					Disposable procedure underpad with virgin fluff and super absorbent polymers, 4 way sealed & waterproof backing. Sizes 60x60cm and 60x90cm		
Unit of Issue		40	35	30	30	30	30	30	30	30	30
Packaging	The packaging has an easily identified description of correct use and size. All necessary detail such as Lot Number and Manufacturer are easy to find	Scores	★ ★ ★ (2.00)		★ ★ ★ (2.00)					★ ★ ★ (2.00)	
	The packaging can be easily opened.	Scores	★ ★ ★ ★ (3.00)		★ ★ ★ ★ (2.80)					★ ★ ★ ★ (2.80)	
Opening and Preparation	The pads are dispensed individually	Scores	★ ★ ★ ★ (2.40)		★ ★ ★ ★ (2.80)					★ ★ ★ ★ (2.40)	
	The packaging remains robust to support the products in a useable condition until empty	Scores	★ ★ ★ ★ (2.80)		★ ★ ★ ★ (2.80)					★ ★ ★ ★ (2.80)	
Clinical Use	The pad in its dry state is comfortable for a procedure	Scores	★ ★ ★ ★ (2.40)		★ ★ ★ ★ (2.00)					★ ★ ★ ★ (2.40)	
	The colour of the absorbent surface does not impede clinical decision making	Scores	★ ★ ★ ★ (2.60)		★ ★ ★ ★ (2.80)					★ ★ ★ ★ (2.80)	
Examples		Sizes	60x60cm	60x90cm	60x40cm	60x60cm	60x75cm	60x90cm	80x180cm	60x60cm	60x90cm
Light absorption required Examinations, insertion of lines, skin prep, dressing & wound procedures, catheter procedures, oral secretions, bottom end procedures, hygiene needs, formed faecal incontinence and endoscopy	Fluid management - light absorbency	Scores	✓	✓	✓	✓	✓	✓	✓	✓	✓
	Fluid management - moderate absorbency	Scores	✗	✓	✗	✓	✓	✓	✓	✓	✓
	Fluid management - high absorbency	Scores	✗	✓	✗	✗	✗	✓	✓	✗	✓
Moderate absorption required Maternity-postpartum, theatre blood loss, wound wash out and oedematous patients / limbs	Surface of the pad to remain dry when pad has been moistened	Scores	★ ★ (0.00) **		★ ★ (0.00) **					★ ★ (0.00) **	
	The pad to remain waterproof to protect the surface underneath it	Scores	★ ★ ★ (2.00) **		★ ★ ★ (2.00) **					★ ★ ★ (1.80) **	
	Retention of liquid within the pad	Scores	★ ★ ★ ★ (1.40)		★ ★ ★ ★ (1.60)					★ ★ ★ ★ (1.80)	
High absorption required Loose faecal incontinence, Maternity-intrapartum, theatre wash-out, theatre trauma blood loss and loss of skin barrier e.g. burns fluid seep		Scores	★ ★ ★ ★ (1.40) **		★ ★ ★ ★ (1.80) **					★ ★ ★ ★ (1.80) **	
	Disposal	Pad can be disposed safely	Scores	★ ★ ★ (1.40) **		★ ★ ★ (1.80) **					★ ★ ★ (1.80) **

**Maximum number of stars attainable 2