




Supplier	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)
Brand	Lang Skintact Tab 'n' Snap	Lang Skintact Tab 'n' Snap	Lang Skintact Tab 'n' Snap	Lang Skintact Tab 'n' Snap	Lang Skintact Tab 'n' Snap
ADULT ECG ELECTRODES					
MPC	FS-WA01C+/10	XL-BM01	FS-WA10C	F-WA10C	FS-VM01/10
NPC	FDK932	FDK912	FDK155	FDK231	FDK420
Minimum Order (UOI Description)	Box	Box	Pouch	Pouch	Box
Quantity of electrodes in minimum order.(UOI)	500	1000	50	50	500
Quantity electrodes in a pouch	10	50	50	50	10
Presentation in pouch S = Single M = Multiple items on backing protection	S	S	M	M	S
Connector Type	Dual Connector Offset	Dual Connector Offset	Dual Connector Onset	Dual Connector Onset	Dual Connector Offset
Gel Type	Wet	Wet	Wet	Wet	Wet
Shape	Teardrop	Teardrop	Rectangle	Teardrop	Teardrop
Size (mm)	37 x 32	45 x 35	37 x 32	37 x 32	45 x 35
X-Ray / MRI	Y	N	Y	Y	N
Shelf Life Open Pouch	7days	7 days	7days	7days	7days
Wear Time	72hours	72 hours	72hours	72hours	72hours
Feature			Enhanced Adhesive		Enhanced Adhesive
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam	FS-WA01-10A (Different Foam) FSWA01C+/50 (50 in Pouch)				FS-VM01 (50 in pouch)
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)
Does the pouch clearly identify product wear time	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★☆☆ (1.90)	★★★☆☆ (1.90)	★★★☆☆ (1.90)	★★★☆☆ (1.90)	★★★☆☆ (1.90)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★☆ (2.30)	★★★★☆ (2.30)	★★★★☆ (2.30)	★★★★☆ (2.30)	★★★★☆ (2.30)
Is the manufacturing product code clearly displayed	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00))
Does the Description clearly identify the electrode properties inside.	★★★☆☆ (1.80)	★★★☆☆ (1.80)	★★★☆☆ (1.80)	★★★☆☆ (1.80)	★★★☆☆ (1.80)
Does the pouch / Pouchet indicate clinical use.	★★★☆☆ (1.20)	★★★☆☆ (1.20)	★★★☆☆ (1.20)	★★★☆☆ (1.20)	★★★☆☆ (1.20)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★☆☆ (1.90)	★★★☆☆ (1.90)	★★★☆☆ (1.90)	★★★☆☆ (1.90)	★★★☆☆ (1.90)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★☆☆ (1.80)	★★★☆☆ (1.80)	★★★☆☆ (1.80)	★★★☆☆ (1.80)	★★★☆☆ (1.80)
Do the electrodes remove easily from the backing piece	★★★☆☆ (2.20)	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)
Does the product apply to the skin correctly and efficiently without wastage.	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)	★★★☆☆ (2.00)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★★☆ (0.00)	★★★★☆ (0.00)	★★★★☆ (0.00)	★★★★☆ (0.00)	★★★★☆ (0.00)
Is the electrode comfortable on the skin for the wearer	★★★☆☆ (1.60)	★★★★☆ (1.80)	★★★☆☆ (2.00)	★★★★☆ (2.20)	★★★☆☆ (2.00)
Any indications / aids for removal on the product	★★★☆☆ (2.00)	★★★★☆ (1.80)	★★★☆☆ (2.00)	★★★★☆ (2.00)	★★★☆☆ (2.00)
Discomfort or skin reaction on removal	★★★☆☆ (1.80)	★★★☆☆ (1.40)	★★★☆☆ (1.60)	★★★☆☆ (2.00)	★★★☆☆ (2.00)
Is the pouch recyclable	★★★★☆ (0.40)	★★★★☆ (0.40)	★★★★☆ (0.40)	★★★★☆ (0.40)	★★★★☆ (0.40)

Supplier	FANNIN (STROUD)	FANNIN (STROUD)	GRAPHIC CONTROLS LTD
Brand	Lang Skintact Tab 'n' Snap	Lang Skintact Tab 'n' Snap	Vermed
ADULT ECG ELECTRODES			
MPC	FS-VM01/A1.5	FS-VF01C	4010429C
NPC	FDK553	FDK342	FJE7450
Minimum Order (UOI Description)	Box	Box	Box
Quantity of electrodes in minimum order.(UOI)	1000	500	1200
Quantity electrodes in a pouch	50	50	60
Presentation in pouch S = Single M = Multiple items on backing protection	S	S	M
Connector Type	Dual Connector Offset	Dual Connector Offset	Dual Connector Onset
Gel Type	Wet	Wet	Wet
Shape	Teardrop	Teardrop	Rectangle
Size (mm)	45 x 35	35 x 45	45 x 32
X-Ray / MRI	N	Y	Y
Shelf Life Open Pouch	7days	7days	30 days
Wear Time	72hours	72hours	48 hours
Feature		Enhanced Adhesive	Resealable Pouch
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam	FSVM01/A1.5/10 (10 in Pouch, 500 in Box)	FS-VF01 (Not XRAY/MRI, minimum order pouch of 50) FS-VF01/10 (10 in Pouch)	
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★ (2.00)	★★★ (2.00)	★★★ (1.80)
Does the pouch clearly identify product wear time	★★★ (2.00)	★★★ (2.00)	★★★ (1.40)
Does the pouch clearly identify the shelf life once pouch opened.	★★★ (1.90)	★★★ (1.90)	★★★ (0.00)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (2.30)	★★★★ (2.30)	★★★ (2.00)
Is the manufacturing product code clearly displayed	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)
Does the Description clearly identify the electrode properties inside.	★★★ (1.80)	★★★ (1.80)	★★★ (1.60)
Does the pouch / Pouchet indicate clinical use.	★★★ (1.20)	★★★ (1.20)	★★★ (0.80)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★ (1.90)	★★★ (1.90)	★★★ (2.00)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★ (1.80)	★★★ (1.80)	★★★ (1.80)
Do the electrodes remove easily from the backing piece	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)
Does the product apply to the skin correctly and efficiently without wastage.	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)
Is the electrode comfortable on the skin for the wearer	★★★★ (2.50)	★★★ (2.00)	★★★ (2.00)
Any indications / aids for removal on the product	★★★ (1.80)	★★★ (2.00)	★★★ (1.60)
Discomfort or skin reaction on removal	★★★ (1.60)	★★★ (1.60)	★★★ (1.60)
Is the pouch recyclable	★★★ (0.40)	★★★ (0.40)	★★★ (1.60)