
















Supplier	ALBERT WAESCHLE LTD	AMBU UK LTD	AMBU UK LTD	CAREFUSION UK 306 LTD	CAREFUSION UK 306 LTD
Brand	Lessa	Ambu White Sensor	Ambu White Sensor	Vital Signs	Silver Trace
ADULT ECG ELECTRODES					
MPC	99.832.02	4535M	4500M	2014777-001	2014786-0001
NPC	FJE6291	FDK047	FDK246	FDK775	FDK777
Minimum Order (UOI Description)	Box	Pouch	Pouch	Box	Box
Quantity of electrodes in minimum order.(UOI)	30	60	30	300	300
Quantity electrodes in a pouch	30	60	30	30	30
Presentation in pouch S = Single M = Multiple items on backing protection	S	M	M	S	M
Connector Type	Onset Stud	Onset Stud	Onset Stud	Onset Stud	Onset Stud
Gel Type	Wet	Wet	Wet	Wet	Wet
Shape	Circular	Rectangle	Teardrop	Circular	Rectangle
Size (mm)	50	44 x 32	50	55	45 x 32
X-Ray / MRI	N	N	Y	N	Y
Shelf Life Open Pouch	7days	1 month	1 month	7 days	7 days
Wear Time	72hours	24hrs	24hrs	72hours	72hours
Feature					
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam		4530M (XRAY/MRI Box 300) 4540M (4mm Banana Adaptor)	4570M (Shape/Size)		2014778-001 (4mm Banana Adaptor - Not XRAY/MRI)
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★★ (1.87)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the pouch clearly identify product wear time	★★★★ (1.50)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★★ (1.60)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (0.00)	★★★★ (0.00)
Is the manufacturing product code clearly displayed	★★★★ (2.00)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (2.00)	★★★★ (2.00)
Does the Description clearly identify the electrode properties inside.	★★★★ (1.20)	★★★★ (1.50)	★★★★ (1.50)	★★★★ (1.40)	★★★★ (1.40)
Does the pouch / Pouchet indicate clinical use.	★★★★ (2.60)	★★★★ (1.00)	★★★★ (1.00)	★★★★ (0.40)	★★★★ (0.40)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★★ (2.00)	★★★★ (1.26)	★★★★ (1.26)	★★★★ (2.00)	★★★★ (2.00)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.60)	★★★★ (1.60)
Do the electrodes remove easily from the backing piece	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the product apply to the skin correctly and efficiently without wastage.	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Is the electrode comfortable on the skin for the wearer	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Any indications / aids for removal on the product	★★★★ (1.00)	★★★★ (1.60)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)
Discomfort or skin reaction on removal	★★★★ (0.20)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.40)	★★★★ (2.00)
Is the pouch recyclable	★★★★ (0.40)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)

Supplier	CAREFUSION UK 306 LTD	ConMed	ConMed	ConMed	Covidien UK
Brand	Silver Trace	Dynatrace	Exactrace	Dynatrace Mini Stress	Kendall
ADULT ECG ELECTRODES					
MPC	900703-805	1500-003	1510-030	1555-005	31.1345.21
NPC	FJE6347		FDK645	FJE6309	FDK714
Minimum Order (UOI Description)	Box	Box	Box	Box	Case
Quantity of electrodes in minimum order.(UOI)	600	600	600	600	300
Quantity electrodes in a pouch	5	3	30	5	30
Presentation in pouch S = Single M = Multiple items on backing protection	M	M	S	M	S
Connector Type	Onset Stud	Onset Stud	Onset Stud	Onset Stud	Onset Stud
Gel Type	Wet	Wet	Wet	Wet	Wet
Shape	Oval	Circular	Circular	Circular	Circular
Size (mm)	50 x 38		51	38	50 x 45
X-Ray / MRI	N	Y	N	Y	N
Shelf Life Open Pouch	7 days				45 days
Wear Time	72hours				72 hours
Feature					
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam	900703-830 (30 Loose in Pouch)				
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★★ (2.00)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (2.00)
Does the pouch clearly identify product wear time	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Does the pouch clearly identify the shelf life once pouch opened.	N/A	N/A	★★★★ (0.00)	N/A	★★★★ (0.00)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Is the manufacturing product code clearly displayed	★★★★ (2.00)	★★★★ (1.73)	★★★★ (1.73)	★★★★ (1.73)	★★★★ (1.65)
Does the Description clearly identify the electrode properties inside.	★★★★ (1.40)	★★★★ (1.30)	★★★★ (1.30)	★★★★ (1.30)	★★★★ (1.70)
Does the pouch / Pouchet indicate clinical use.	★★★★ (0.40)	★★★★ (1.20)	★★★★ (1.20)	★★★★ (1.20)	★★★★ (1.00)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	N/A	N/A	★★★★ (1.80)	N/A	★★★★ (1.80)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	N/A	N/A	★★★★ (1.60)	N/A	★★★★ (1.60)
Do the electrodes remove easily from the backing piece	★★★★ (1.40)	★★★★ (2.00)	★★★★ (1.80)	★★★★ (2.00)	★★★★ (1.40)
Does the product apply to the skin correctly and efficiently without wastage.	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.80)	★★★★ (2.00)	★★★★ (1.80)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★★ (0.00)	★★★★ (1.40)	★★★★ (1.20)	★★★★ (1.40)	★★★★ (0.00)
Is the electrode comfortable on the skin for the wearer	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.20)	★★★★ (1.80)
Any indications / aids for removal on the product	★★★★ (0.60)	★★★★ (0.40)	★★★★ (1.40)	★★★★ (0.00)	★★★★ (1.80)
Discomfort or skin reaction on removal	★★★★ (1.40)	★★★★ (2.20)	★★★★ (1.40)	★★★★ (1.60)	★★★★ (1.40)
Is the pouch recyclable	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.50)

Supplier	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)
Brand	Lang Skintact Premier	Lang Skintact Premier	Lang Skintact Premier	Lang Skintact Premier	Lang Skintact Premier
ADULT ECG ELECTRODES					
MPC	FS-TF	FS-RE/6	FS-50C	FS-RB5	FS-TA
NPC	FDK234	FDK121	FDK341	FDK116	FDK233
Minimum Order (UOI Description)	Pouch	Pouch	Box	Pouch	Pouch
Quantity of electrodes in minimum order.(UOI)	30	30	300	30	30
Quantity electrodes in a pouch	30	30	30	30	30
Presentation in pouch S = Single M = Multiple items on backing protection	M	M	S	M	S
Connector Type	Onset Stud	Onset Stud	Onset Stud	4mm banana adaptor	Onset Stud
Gel Type	Wet	Wet	Wet	Wet	Wet
Shape	Teardrop	Rectangle	Circular	Rectangle	Circular
Size (mm)	45	44 x 32	50	44 x 28	37
X-Ray / MRI	N	N	Y	N	N
Shelf Life Open Pouch	7 days	7 days	7 days	7 days	7 days
Wear Time	72hours	72 hours	72hours	72 hours	72 hours
Feature					Enhanced Adhesive
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam	FS-TF/6 (6 electrodes to backing strip in a case of 1200)				
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the pouch clearly identify product wear time	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★★ (0.00)	★★★★ (1.00)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (2.50)	★★★★ (2.50)	★★★★ (2.30)	★★★★ (2.30)	★★★★ (2.30)
Is the manufacturing product code clearly displayed	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the Description clearly identify the electrode properties inside.	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)
Does the pouch / Pouchet indicate clinical use.	★★★★ (1.20)	★★★★ (1.20)	★★★★ (1.20)	★★★★ (1.20)	★★★★ (1.20)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)
Do the electrodes remove easily from the backing piece	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.83)
Does the product apply to the skin correctly and efficiently without wastage.	★★★★ (2.00)	★★★★ (2.50)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Is the electrode comfortable on the skin for the wearer	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.80)	★★★★ (1.60)	★★★★ (1.83)
Any indications / aids for removal on the product	★★★★ (2.00)	★★★★ (2.30)	★★★★ (1.80)	★★★★ (2.00)	★★★★ (1.83)
Discomfort or skin reaction on removal	★★★★ (1.83)	★★★★ (1.80)	★★★★ (1.20)	★★★★ (1.80)	★★★★ (1.83)
Is the pouch recyclable	★★★★ (0.40)	★★★★ (0.40)	★★★★ (0.40)	★★★★ (0.40)	★★★★ (0.40)

Supplier	FANNIN (STROUD)	FANNIN (STROUD)	GRAPHIC CONTROLS LTD	MEDIMAXTECH UK LTD	PHILIPS HEALTHCARE
Brand	Lang Skintact Premier	Lang Skintact Premier	Vermed	MMT	Philips Healthcare
ADULT ECG ELECTRODES					
MPC	F-TB	F-REC	FW50	MXF 55RW	M2202A
NPC	FDK124	FDK123	FJE6299	FJE6339	FDI237
Minimum Order (UOI Description)	Pouch	Pouch	Box	Case	Box
Quantity of electrodes in minimum order.(UOI)	30	30	400	6000	300
Quantity electrodes in a pouch	30	30	50	30	5
Presentation in pouch S = Single M = Multiple items on backing protection	S	M	S	S	M
Connector Type	Onset Stud	Onset Stud	Onset Stud	Onset Stud	Onset Stud
Gel Type	Wet	Wet	Wet	Wet	Wet
Shape	Circular	Rectangle	Circular	Circular	Circular
Size (mm)	50	44 x 32	50	55	50
X-Ray / MRI	N	Y	N	N	Y
Shelf Life Open Pouch	7 days	7days	15 days		7 days
Wear Time	72hours	72hours	72 hours		3 days
Feature					Enhanced Adhesive
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam	FTB/3 Quantity in Pouch		FW50R (XRAY/MRI)		
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.60)	★★★★ (2.00)
Does the pouch clearly identify product wear time	★★★★ (2.00)	★★★★ (2.00)	★★★ (0.50)	★★★ (0.50)	★★★ (0.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★★ (1.90)	★★★★ (1.90)	★★★ (0.00)	★★★ (0.00)	N/A
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (2.30)	★★★★ (2.30)	★★★★ (1.90)	★★★★ (1.20)	★★★★ (2.30)
Is the manufacturing product code clearly displayed	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.10)	★★★★ (1.80)
Does the Description clearly identify the electrode properties inside.	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.60)	★★★★ (1.00)
Does the pouch / Pouchet indicate clinical use.	★★★★ (1.20)	★★★★ (1.20)	★★★★ (2.40)	★★★★ (1.70)	★★★★ (1.00)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★★ (1.90)	★★★★ (1.90)	★★★ (0.30)	★★★ (0.80)	N/A
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.60)	★★★★ (1.70)	N/A
Do the electrodes remove easily from the backing piece	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.60)	★★★★ (1.60)	★★★★ (2.00)
Does the product apply to the skin correctly and efficiently without wastage.	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.25)	★★★★ (2.00)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)
Is the electrode comfortable on the skin for the wearer	★★★★ (2.00)	★★★★ (1.80)	★★★★ (2.00)	★★★★ (1.60)	★★★★ (2.00)
Any indications / aids for removal on the product	★★★★ (2.00)	★★★★ (1.80)	★★★ (0.00)	★★★★ (1.80)	★★★★ (2.00)
sDiscomfort or skin reaction on removal	★★★★ (1.83)	★★★★ (1.80)	★★★★ (1.40)	★★★★ (1.60)	★★★★ (1.80)
Is the pouch recyclable	★★★ (0.40)	★★★ (0.40)	★★★ (0.00)	★★★ (0.40)	★★★ (0.00)

Supplier	PHILIPS HEALTHCARE	PHILIPS HEALTHCARE
Brand	Philips Healthcare	Philips Healthcare
ADULT ECG ELECTRODES		
MPC	40420A	40493E
NPC	FJE6380	FDK895
Minimum Order (UOI Description)	Box	Box
Quantity of electrodes in minimum order.(UOI)	1000	300
Quantity electrodes in a pouch	30	30
Presentation in pouch S = Single M = Multiple items on backing protection	M	S
Connector Type	Onset Stud	Onset Stud
Gel Type	Wet	Wet
Shape	Teardrop	Circular
Size (mm)	65 x 28	54
X-Ray / MRI	N	N
Shelf Life Open Pouch	7 days	7 days
Wear Time	3 days	3 days
Feature	Enhanced Adhesive	
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam		40493D (5 in pouch)
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★ (2.00)	★★★ (2.00)
Does the pouch clearly identify product wear time	★★★ (0.00)	★★★ (0.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★ (0.00)	★★★ (0.00)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★ (1.60)	★★★ (1.60)
Is the manufacturing product code clearly displayed	★★★ (1.80)	★★★ (1.80)
Does the Description clearly identify the electrode properties inside.	★★★ (1.40)	★★★ (1.40)
Does the pouch / Pouchet indicate clinical use.	★★★ (1.00)	★★★ (1.00)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★ (2.00)	★★★ (2.00)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★ (1.80)	★★★ (1.80)
Do the electrodes remove easily from the backing piece	★★★ (1.80)	★★★ (2.00)
Does the product apply to the skin correctly and efficiently without wastage.	★★★ (2.00)	★★★ (2.00)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★ (0.00)	★★★ (0.00)
Is the electrode comfortable on the skin for the wearer	★★★ (2.00)	★★★ (2.00)
Any indications / aids for removal on the product	★★★ (1.80)	★★★ (1.80)
Discomfort or skin reaction on removal	★★★ (2.00)	★★★ (1.80)
Is the pouch recyclable	★★★ (0.00)	★★★ (0.00)