





















Supplier	3M UK PLC	ALBERT WAESCHLE LTD	AMBU UK LTD	AMBU UK LTD	AMBU UK LTD
Brand	3M Red Dot	Lessa	Ambu Blue Sensor	Ambu White Sensor	Ambu Blue Sensor
ADULT ECG ELECTRODES					
MPC	2330	99.832.38	SU-00-C/100	0715M	BS-Tab-2300
NPC	FDI536	FJE6292	FDK061	FDI307	FDI306
Minimum Order (UOI Description)	Pouch	Box	Pouch	Box	Box
Quantity of electrodes in minimum order.(UOI)	100	100	100	500	500
Quantity electrodes in a pouch	100	100	100	100	100
Presentation in pouch S = Single M = Multiple items on backing protection	M	M	M	M	M
Connector Type	Tab	Tab	Tab	Tab	Tab
Gel Type	Solid	Solid	Wet	Solid	Solid
Shape	Rectangle	Rectangle	Teardrop	Rectangle	Rectangle
Size (mm)	20 x 30	34 x 22	49 x 33	34 x 26	34 x 23
X-Ray / MRI	Y	N	N	Y	Y
Shelf Life Open Pouch	7 days		1 month	1 month	1 month
Wear Time			24hrs	2 hrs	2 hrs
Feature			Enhanced Adhesive		Enhanced Adhesive
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam			SU-00-A (Banana Adaptor)	0615M (size 37x26) 0415M (5000 in box, size 34x23) 0315M (size 35x21)	
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★★ (1.80)	★★★★ (1.87)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the pouch clearly identify product wear time	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (0.00)	★★★★ (1.00)	★★★★ (2.00)	★★★★ (2.20)	★★★★ (2.20)
Is the manufacturing product code clearly displayed	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.90)	★★★★ (1.80)	★★★★ (1.80)
Does the Description clearly identify the electrode properties inside.	★★★★ (2.00)	★★★★ (1.20)	★★★★ (1.50)	★★★★ (1.60)	★★★★ (1.60)
Does the pouch / Pouchet indicate clinical use.	★★★★ (1.50)	★★★★ (1.60)	★★★★ (1.00)	★★★★ (1.00)	★★★★ (1.00)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.26)	★★★★ (1.80)	★★★★ (1.80)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★★ (1.60)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.20)	★★★★ (1.20)
Do the electrodes remove easily from the backing piece	★★★★ (2.00)	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.30)	★★★★ (2.30)
Does the product apply to the skin correctly and efficiently without wastage.	★★★★ (2.20)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★★ (1.20)	★★★★ (1.20)	★★★★ (2.20)	★★★★ (1.00)	★★★★ (1.00)
Is the electrode comfortable on the skin for the wearer	★★★★ (2.20)	★★★★ (2.20)	★★★★ (2.00)	★★★★ (2.20)	★★★★ (2.40)
Any indications / aids for removal on the product	★★★★ (1.60)	★★★★ (1.60)	★★★★ (1.40)	★★★★ (1.60)	★★★★ (2.20)
Discomfort or skin reaction on removal	★★★★ (2.20)	★★★★ (1.20)	★★★★ (1.80)	★★★★ (2.20)	★★★★ (2.60)
Is the pouch recyclable	★★★★ (0.50)	★★★★ (0.40)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)

Supplier	CAREFUSION UK 306 LTD	CAREFUSION UK 306 LTD	ConMed	ConMed	Covidien UK
Brand	Vital Signs	Vital Signs	Fastrace	HeartBeat	Kendall
ADULT ECG ELECTRODES					
MPC	2009101-003P	9623-810P	1915-100	1925-100	31433538
NPC	FDK1005	FDK941	FDK651	FDK653	FDK877
Minimum Order (UOI Description)	Box	Box	Box	Box	Case
Quantity of electrodes in minimum order.(UOI)	1000	1000	1000	1000	4000
Quantity electrodes in a pouch	100	100	100	100	100
Presentation in pouch S = Single M = Multiple items on backing protection	M	M	M	M	M
Connector Type	Tab	Tab	Tab	Tab	Tab
Gel Type	Solid	Solid	Solid	Solid	Solid
Shape	Rectangle	Rectangle	Rectangle	Rectangle	Rectangle
Size (mm)	34 x 26	34 x 22	34 x 22	32 x 20	33 x 23
X-Ray / MRI	N	N	N	N	N
Shelf Life Open Pouch	7 days				45 days
Wear Time	1 hour	1 hour			72 hours
Feature					
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam	2009101-003 (adhesive change, size 34x23)				
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (2.00)
Does the pouch clearly identify product wear time	★★★ (0.00)	★★★★ (1.00)	★★★ (0.00)	★★★ (0.00)	★★★★ (1.50)
Does the pouch clearly identify the shelf life once pouch opened.	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★ (0.00)	★★★★ (2.00)	★★★ (0.00)	★★★ (0.00)	★★★★ (2.30)
Is the manufacturing product code clearly displayed	★★★★ (1.80)	★★★★ (2.00)	★★★★ (1.73)	★★★★ (1.73)	★★★★ (1.80)
Does the Description clearly identify the electrode properties inside.	★★★★ (1.25)	★★★★ (1.25)	★★★★ (1.30)	★★★★ (1.30)	★★★★ (1.40)
Does the pouch / Pouchet indicate clinical use.	★★★★ (1.25)	★★★★ (1.25)	★★★★ (1.20)	★★★★ (1.20)	★★★★ (2.00)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★★ (1.00)	★★★★ (2.00)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★★ (1.80)	★★★★ (1.60)	★★★★ (1.60)	★★★★ (1.60)	★★★★ (1.70)
Do the electrodes remove easily from the backing piece	★★★★ (2.00)	★★★★ (1.80)	★★★★ (1.40)	★★★★ (2.50)	★★★★ (1.80)
Does the product apply to the skin correctly and efficiently without wastage.	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.40)	★★★★ (2.00)	★★★★ (2.00)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★ (0.00)	★★★★ (1.20)	★★★★ (2.60)	★★★★ (2.40)	★★★ (1.20)
Is the electrode comfortable on the skin for the wearer	★★★★ (2.00)	★★★★ (2.20)	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.00)
Any indications / aids for removal on the product	★★★★ (2.00)	★★★★ (2.00)	★★★ (1.20)	★★★★ (1.60)	★★★★ (2.00)
Discomfort or skin reaction on removal	★★★★ (2.00)	★★★★ (2.00)	★★★ (1.20)	★★★★ (2.00)	★★★★ (1.80)
Is the pouch recyclable	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)	★★★ (0.50)

Supplier	Covidien UK	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)
Brand	Kendall	Lang Skintact EasiTabs	Lang Skintact EasiTabs	Lang Skintact EasiTabs	Lang Skintact EasiTabs
ADULT ECG ELECTRODES					
MPC	30807732	FS-WB00	F-WA00	RT-38	RT74BD
NPC	FDK731	FDK244	FDK117	FDK245	FDK1053
Minimum Order (UOI Description)	Case	Box	Pouch	Box	Case
Quantity of electrodes in minimum order.(UOI)	2000	500	50	500	5000
Quantity electrodes in a pouch	100	100	50	100	100
Presentation in pouch S = Single M = Multiple items on backing protection	M	M	M	M	M
Connector Type	Tab	Tab	Tab	Tab	Tab
Gel Type	Solid	Wet	Wet	Solid	Solid
Shape	Rectangle	Square	Rectangle	Rectangle	Rectangle
Size (mm)	33 x 23	28	37 x 32	25 x 21	26 x 24
X-Ray / MRI	N	Y	Y	N	N
Shelf Life Open Pouch	45 days	7 days	7 days	7 days	7 days
Wear Time	72 hours	72hours	72hours	24 hours	24 hours
Feature					
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam					RT74 (Minimum Order of 1 Pouch)
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)
Does the pouch clearly identify product wear time	★★★ (1.50)	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)	★★★★ (2.40)	★★★★ (2.40)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (2.30)	★★★★ (2.50)	★★★★ (2.50)	★★★★ (2.50)	★★★★ (2.50)
Is the manufacturing product code clearly displayed	★★★ (1.80)	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)
Does the Description clearly identify the electrode properties inside.	★★★ (1.40)	★★★ (1.80)	★★★ (1.80)	★★★ (1.80)	★★★ (1.80)
Does the pouch / Pouchet indicate clinical use.	★★★ (2.00)	★★★ (1.20)	★★★ (1.20)	★★★ (1.50)	★★★ (1.50)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★ (1.80)	★★★ (1.80)	★★★ (1.80)	★★★ (1.80)	★★★ (1.80)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★ (1.70)	★★★ (1.80)	★★★ (1.80)	★★★ (1.80)	★★★ (1.80)
Do the electrodes remove easily from the backing piece	★★★ (1.80)	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)
Does the product apply to the skin correctly and efficiently without wastage.	★★★ (2.00)	★★★ (1.80)	★★★ (2.00)	★★★ (2.00)	★★★ (2.00)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★ (1.20)	★★★ (1.80)	★★★ (1.80)	★★★ (1.20)	★★★ (1.20)
Is the electrode comfortable on the skin for the wearer	★★★ (2.20)	★★★ (2.20)	★★★ (2.20)	★★★★ (2.40)	★★★ (2.20)
Any indications / aids for removal on the product	★★★ (1.80)	★★★ (2.20)	★★★ (2.20)	★★★ (2.00)	★★★ (1.80)
Discomfort or skin reaction on removal	★★★ (2.00)	★★★★ (2.40)	★★★ (2.20)	★★★ (1.80)	★★★ (1.80)
Is the pouch recyclable	★★★ (0.50)	★★★ (0.40)	★★★ (0.40)	★★★ (0.40)	★★★ (0.40)

Supplier	GRAPHIC CONTROLS LTD	GRAPHIC CONTROLS LTD	MEDIMAXTECH UK LTD	MEDIMAXTECH UK LTD	PHILIPS HEALTHCARE
Brand	Vermed	Vermed	MMT	MMT	Philips Healthcare
ADULT ECG ELECTRODES					
MPC	A10201	4010423C	MXT-2034	MXT-2634	13943B
NPC	FJE6323	FJE7453	FDK983	FDK985	FDK534
Minimum Order (UOI Description)	Box	Box	Case	Case	Box
Quantity of electrodes in minimum order.(UOI)	100	1500	5000	5000	1000
Quantity electrodes in a pouch	100	100	100	100	100
Presentation in pouch S = Single M = Multiple items on backing protection	M	M	M	M	M
Connector Type	Tab	Tab	Tab	Tab	Tab
Gel Type	Solid	Wet	Solid	Solid	Solid
Shape	Rectangle	Teardrop	Rectangle	Rectangle	Rectangle
Size (mm)	34 x 23	52 x 35	34 x 20	34 x 26	35 x 22
X-Ray / MRI	Y	Y	N	N	N
Shelf Life Open Pouch	15 days	30 days			7 days
Wear Time	24 hours	48 hours			3 days
Feature		Resealable Pouch			
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam					
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★★ (1.40)	★★★★ (1.80)	★★★★ (1.40)	★★★★ (1.40)	★★★★ (2.00)
Does the pouch clearly identify product wear time	★★★★ (1.10)	★★★★ (1.40)	★★★★ (1.00)	★★★★ (1.00)	★★★★ (0.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (2.60)	★★★★ (2.00)	★★★★ (1.20)	★★★★ (1.20)	★★★★ (0.00)
Is the manufacturing product code clearly displayed	★★★★ (2.20)	★★★★ (2.00)	★★★★ (1.10)	★★★★ (1.10)	★★★★ (1.80)
Does the Description clearly identify the electrode properties inside.	★★★★ (2.00)	★★★★ (1.60)	★★★★ (1.40)	★★★★ (1.40)	★★★★ (1.40)
Does the pouch / Pouchet indicate clinical use.	★★★★ (2.40)	★★★★ (0.80)	★★★★ (2.40)	★★★★ (2.40)	★★★★ (1.00)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.60)	★★★★ (1.60)	★★★★ (2.00)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)
Do the electrodes remove easily from the backing piece	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)
Does the product apply to the skin correctly and efficiently without wastage.	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.60)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★★ (1.60)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (1.00)
Is the electrode comfortable on the skin for the wearer	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Any indications / aids for removal on the product	★★★★ (1.80)	★★★★ (1.60)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (2.00)
Discomfort or skin reaction on removal	★★★★ (2.00)	★★★★ (1.60)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Is the pouch recyclable	★★★★ (0.00)	★★★★ (1.60)	★★★★ (0.40)	★★★★ (0.40)	★★★★ (0.00)