













Supplier	3M UK PLC	AMBU UK LTD	AMBU UK LTD	AMBU UK LTD	AMBU UK LTD
Brand	3M Red Dot	Ambu Blue Sensor	Ambu Blue Sensor	Ambu Blue Sensor	Ambu White Sensor
ADULT ECG ELECTRODES					
MPC	2239	VLC-00-s/10	VL-00-S	L-00-S	7841P
NPC	FDK148	FDK432	FDK587	FDK188	FDI298
Minimum Order (UOI Description)	Pouch	Pouch	Pouch	Pouch	Pouch
Quantity of electrodes in minimum order.(UOI)	50	10	25	25	60
Quantity electrodes in a pouch	50	10	25	25	60
Presentation in pouch S = Single M = Multiple items on backing protection	S	S	S	S	M
Connector Type	Onset Stud	Offset Stud	Offset Stud	Offset Stud	Onset Stud
Gel Type	Solid	Wet	Wet	Wet	Solid
Shape	Circular	Teardrop	Teardrop	Teardrop	Square
Size (mm)	61	72 x 68	72 x 68	68.2 x 55	42
X-Ray / MRI	N	N	N	N	Y
Shelf Life Open Pouch	30 days	1 month	1 month	1 month	1 month
Wear Time	3 days	72hrs	72hrs	72hrs	24hrs
Feature	Removal Tab	Enhanced Adhesive	Enhanced Adhesive	Enhanced Adhesive	Clear Electrode
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam		VLC00-S/25 (25 in Pouch)			
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the pouch clearly identify product wear time	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★★ (2.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (0.80)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Is the manufacturing product code clearly displayed	★★★★ (2.00)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)
Does the Description clearly identify the electrode properties inside.	★★★★ (2.00)	★★★★ (1.50)	★★★★ (1.50)	★★★★ (1.50)	★★★★ (1.50)
Does the pouch / Pouchet indicate clinical use.	★★★★ (1.50)	★★★★ (1.00)	★★★★ (1.00)	★★★★ (1.00)	★★★★ (1.00)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★★ (2.00)	★★★★ (1.26)	★★★★ (1.26)	★★★★ (1.26)	★★★★ (1.26)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)
Do the electrodes remove easily from the backing piece	★★★★ (2.00)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.80)	★★★★ (2.25)
Does the product apply to the skin correctly and efficiently without wastage.	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.25)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★★ (1.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (2.20)
Is the electrode comfortable on the skin for the wearer	★★★★ (1.80)	★★★★ (1.40)	★★★★ (1.40)	★★★★ (1.60)	★★★★ (2.75)
Any indications / aids for removal on the product	★★★★ (2.00)	★★★★ (1.60)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (2.00)
Discomfort or skin reaction on removal	★★★★ (1.80)	★★★★ (1.00)	★★★★ (1.20)	★★★★ (1.40)	★★★★ (2.30)
Is the pouch recyclable	★★★★ (0.50)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)

Supplier	ConMed	ConMed	Covidien UK	FANNIN (STROUD)	FANNIN (STROUD)
Brand	Cleartrace	Cleartrace LT	Kendall	Lang Skintact Premier	Lang Skintact XL T-V002
ADULT ECG ELECTRODES					
MPC	1700-030	2710-030	22700	XL-TVO01	XL-TVO02
NPC	FDK640	FDK657	FDK707	FDK905	FDK911
Minimum Order (UOI Description)	Box	Box	Case	Case	Case
Quantity of electrodes in minimum order.(UOI)	600	600	1000	300	300
Quantity electrodes in a pouch	30	30	50	30	30
Presentation in pouch S = Single M = Multiple items on backing protection	S	S	M	S	S
Connector Type	Onset Stud	Onset Stud	Onset Stud	Offset Stud	Offset Stud
Gel Type	Solid	Solid	Solid	Wet	Wet
Shape	Circular	Square	Circular	Circular	Circular
Size (mm)	45	50	45	60	68
X-Ray / MRI	N	N	N	N	N
Shelf Life Open Pouch			45 days	7days	7days
Wear Time		7 Days	72hours	72hours	72hours
Feature					
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam		2700-030 (XRAY/MRI)	22750 (XRAY/MRI)		
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★★ (1.90)	★★★★ (1.90)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the pouch clearly identify product wear time	★★★★ (0.00)	★★★★ (0.00)	★★★★ (1.50)	★★★★ (2.00)	★★★★ (2.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.00)	★★★★ (1.90)	★★★★ (1.90)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (0.00)	★★★★ (0.00)	★★★★ (2.30)	★★★★ (2.30)	★★★★ (2.30)
Is the manufacturing product code clearly displayed	★★★★ (1.73)	★★★★ (1.73)	★★★★ (1.80)	★★★★ (2.00)	★★★★ (2.00)
Does the Description clearly identify the electrode properties inside.	★★★★ (1.30)	★★★★ (1.30)	★★★★ (1.60)	★★★★ (1.80)	★★★★ (1.80)
Does the pouch / Pouchet indicate clinical use.	★★★★ (1.20)	★★★★ (1.20)	★★★★ (1.00)	★★★★ (1.20)	★★★★ (1.20)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.90)	★★★★ (1.90)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★★ (1.60)	★★★★ (1.60)	★★★★ (1.60)	★★★★ (1.80)	★★★★ (1.80)
Do the electrodes remove easily from the backing piece	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the product apply to the skin correctly and efficiently without wastage.	★★★★ (2.00)	★★★★ (2.20)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.80)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★★ (0.60)	★★★★ (0.50)	★★★★ (1.60)	★★★★ (0.00)	★★★★ (0.00)
Is the electrode comfortable on the skin for the wearer	★★★★ (2.20)	★★★★ (2.20)	★★★★ (2.00)	★★★★ (1.20)	★★★★ (1.20)
Any indications / aids for removal on the product	★★★★ (1.80)	★★★★ (2.00)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)
Discomfort or skin reaction on removal	★★★★ (1.80)	★★★★ (1.80)	★★★★ (2.00)	★★★★ (1.20)	★★★★ (1.00)
Is the pouch recyclable	★★★★ (0.00)	★★★★ (0.00)	★★★★ (0.50)	★★★★ (0.40)	★★★★ (0.40)

Supplier	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)	FANNIN (STROUD)	PHILIPS HEALTHCARE
Brand	Lang Skintact Elite	Lang Skintact Elite	Lang Skintact Premier	Lang Skintact Premier	Lang Skintact Premier	Philips Healthcare
ADULT ECG ELECTRODES						
MPC	T-VO02	T-VO01	T-602	T-60	CT-601	40489E
NPC	FDK929	FDK902	FDK105	FDK115	FDK108	FDK1003
Minimum Order (UOI Description)	Box	Box	Pouch	Pouch	Pouch	Box
Quantity of electrodes in minimum order.(UOI)	300	300	30	30	30	300
Quantity electrodes in a pouch	30	30	30	30	30	30
Presentation in pouch S = Single M = Multiple items on backing protection	S	S	S	S	S	S
Connector Type	Offset Stud	Offset Stud	Onset Stud	Onset Stud	Onset Stud	Onset Stud
Gel Type	Wet	Wet	Tac	Wet	Tac	Wet
Shape	Circular	Circular	Circular	Circular	Circular	Circular
Size (mm)	72	66 x 60	60	60	50	59
X-Ray / MRI	N	N	N	N	N	N
Shelf Life Open Pouch	7 days	7 days	7 days	7 days	7days	7 days
Wear Time	72hours	72hours	72hours	72hours	72hours	3 days
Feature					Clear Tape Backing	
Variations e.g. Radiotranslucent Enhanced Adhesive Different Foam						
Is the Expiry Date & Lot Number or Manufacturing Date & Shelf Life Clearly visible	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the pouch clearly identify product wear time	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★ (0.00)
Does the pouch clearly identify the shelf life once pouch opened.	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)	★★★ (0.00)
Are there clear instructions for use, including skin-preparation and electrode preparation	★★★★ (2.30)	★★★★ (2.30)	★★★★ (2.30)	★★★★ (2.30)	★★★★ (2.30)	★★★★ (2.30)
Is the manufacturing product code clearly displayed	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (1.80)
Does the Description clearly identify the electrode properties inside.	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.40)
Does the pouch / Pouchet indicate clinical use.	★★★ (1.20)	★★★ (1.20)	★★★ (1.20)	★★★ (1.20)	★★★ (1.20)	★★★ (1.00)
Does the pouch / Pouchet have a space to record date it was opened or new expiry date	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (1.90)	★★★★ (2.00)
Can the pouch / Pouchet be easily opened without splitting or damaging the Packaging using hands only	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)	★★★★ (1.80)
Do the electrodes remove easily from the backing piece	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the product apply to the skin correctly and efficiently without wastage.	★★★★ (2.00)	★★★★ (1.83)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)	★★★★ (2.00)
Does the backing card have a colour or pattern which would make it visible on a floor or in the bedsheets.	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)	★★★ (0.00)
Is the electrode comfortable on the skin for the wearer	★★★ (1.20)	★★★★ (1.33)	★★★★ (1.40)	★★★★ (1.40)	★★★★ (1.40)	★★★★ (1.80)
Any indications / aids for removal on the product	★★★★ (1.60)	★★★★ (1.67)	★★★★ (1.80)	★★★★ (1.60)	★★★★ (1.80)	★★★★ (1.80)
Discomfort or skin reaction on removal	★★★ (1.00)	★★★ (1.17)	★★★★ (1.40)	★★★★ (1.40)	★★★ (1.20)	★★★★ (1.60)
Is the pouch recyclable	★★★ (0.40)	★★★ (0.40)	★★★ (0.40)	★★★ (0.40)	★★★ (0.40)	★★★ (0.00)